

# **MOVE**

*It's your move!*



**MOTION**  
COMPOSITES



# MOVE



## Now it's your move!

When we sat down to create the Move, we began with our usual accept-no-compromises goal — to build the lightest, most efficient wheelchair in its category. Period. It offers the design ingenuity you'd expect from Motion Composites. And it adds something new — accessibility. We are proud to bring Motion Composites superb quality to the standard wheelchair category.

## Simple. Light. Tough.

We designed the Move for lightness and simplicity. Innovative engineering and lightweight accessories have enabled us to create the lightest standard wheelchair on the market. And, because we want to maintain a lightweight package, we offer Newton One spoke wheels standard on every Move.

## Technical Specifications

### Structure

Frame	Folding with aluminum crossbrace
Material	Aluminum alloy
Weight	16 lbs. (7.25 kg) (w/o rear wheels & footrest)
Weight capacity	250 lbs. (113 kg) Kit HD: 350 lbs. (159 kg)

### Dimensions

Width even sizes	14" (35.6 cm) to 22" (55.9 cm) in
Depth (50.8 cm)	16" (40.6 cm), 18" (45.7 cm), 20"
Front seat-to-floor height	12 ½" (31.8 cm) to 20 ½" (52.1 cm)
Rear seat-to-floor height	12 ¼" (31.1 cm) to 20" (50.8 cm)

70294-a.0-MOVE FLYER - ENGLISH