PRODUCT



LIGHTWEIGHT & COMPACT •

Upright



• 10° Fixed Tilt

Provides upright posture to enhance functional positioning

- Self-Tensio™

 Built in self-tensioning seating system increases sitting tolerance
- 5 Years of Growth Capabilities
 Extends chair life for a growing child
- Comprehensive Positioning Options
 Provides maximum comfort and optimal positioning

Models











Vivo

EZ Rider

EZ Planar

Convertible

CV Planar

SAFE & DEPENDABLE •

Fixed Tilt



- 30° Fixed Tilt
 - Increased tilt angle improves postural control through upper body positioning
- **Self-Tensio™**Built in self-tensioning seating system increases sitting tolerance
- 5 Years of Growth Capabilities
 Extends chair life for a growing child
- Comprehensive Positioning Options
 Provides maximum comfort and optimal positioning

Models









Metro

Cruiser

CX Planar

Scout

CUSTOMIZED & CONVENIENT

Tilt-in-Space



• Up to 45° Adjustable Tilt

Provides pressure relief and upper body support

Dynamic Positioning

Accommodates functional or rest positions and increased sitting tolerance for maximum comfort

• Adjustable Recline

Opens the seat-to-back angle to provide postural relief

• 5 Years of Growth Capabilities

Extends chair life for a growing child









CuddleBug

Rodeo

Safari

Child Restraint System



Carrot 3

- Adaptive Positioning
 High density foam inserts
 for optimal positioning and
 maximum comfort
- Grows With Your Child 30-108 lbs / 37" 60"
- "Free Angle" Recline 180
 Allows for easy diaper changes and conforms to vehicle's seat
- LATCH Equipped Lower Anchors and Tethers for Children
- Aircraft Travel Certified
 Making travel easier for the whole family

Accessories









BEST IN CLASS

Convaid mobility products are engineered for superior performance and safety. Our rigorous internal testing far exceeds industry standards and FDA requirements. Every product is backed by an unconditional lifetime warranty.

CUSTOMIZED

Each wheelchair is individually assembled to order. We partner with dealers, therapists and caregivers to customize each wheelchair to meet one young person's individual physical needs.

COMPASSIONATE

Convaid is dedicated to helping special needs families reclaim their freedom. We understand that families and caregivers need a wheelchair that is lightweight, convenient and easy to handle, colorful and contemporary in style.



Now You're Going Places



888-Convaid (266-8243) www.convaid.com

