

CATALYST 5 & 5T





INNOVATIVE AND PATENTED FRAME DESIGN

Patented frame concept uses intelligent geometry to create a significantly stronger, lighter and better performing chair.

STANDARD & CURVED AXLE PLATES

Catalyst 5 provides the ability to adjust CG over a 3 1/2" range. The adjustment is a snap with the innovative 3D locking index washer design and requires only one wrench.

SWING AWAY FOOTRESTS

Enjoy the simplicity and strength of Catalyst's patented swing-away footrest hangers. The tapered engagement keeps the footrests snug and rattle free. Comes with a lifetime warranty.

UNMATCHED PERFORMANCE.

The professional's best choice in folding ultra-lightweight mobility.

Also Available in Titanium!



SPECIFICATIONS

5 Transport Weight:	12.9 lbs.*	5.85 kg*
5Ti Transport Weight:	10.9 lbs.*	4.94 kg*
Seat Width:	14-22 in.	36-56 cm
Seat Depth:	14-20 in.	36-51 cm
Front Seat Height:	13.5-21 in.	34-53 cm
Rear Seat Height:	11.5-19.5 in.	29-50 cm
Top Frame Tubing Size:	1 in.	2.54 cm
Bottom Frame Tubing Size:	1.125 in.	2.86 cm
Weight Capacity:	300 lbs.	136 kg
Ansi Resna WC-19 Transit:	Yes	Yes

* 16"x16" Frame, 5"x1" Poly Casters, Swing Away Arm Receivers. Footrest and Rear Wheels removed.