

CATALYST X



5VX

IGNITING THE EVOLUTION

ENDLESS **STYLE**. EXPECTED **QUALITY**. ESSENTIAL **VALUE**.



From the logo to the frame, Catalyst 5VX has been redesigned with a modern, minimalist approach. Enhanced designs let you optimize your ride.

THE RIDE YOU'VE COME TO EXPECT FROM CATALYST.

BETTER STYLE


Be bold. Be you. With thousands of combinations between frame color, anodizing package, options and accessories, your style begins with you.



BETTER ADJUSTMENT

Center of Gravity adjustment and ride performance go hand in hand. The redesigned axle plate and rear frame give you more adjustment over a greater range. Find your fit, ease your propulsion and optimize the ride you've come to expect with Catalyst.

Adjustable Camber gives you up to 4° of camber in 1" increments. A simple adjustment and convenient option you can optimize now or later without any additional parts.



MORE OPTIONS & ACCESSORIES

SWING-IN ANTI-TIPS

The new swing in anti-tips make it easy for the user to operate. By simply pushing down on the ergonomically designed handle you can swing them in, out of the way.



TRANSIT WHEELS

Integrated into the rear frame, Transit Wheels allow you to remove your rear wheels for easier access through narrow spaces.



FENDER SIDE GUARDS

Choose between carbon fiber and foam padded aluminum side guards with fender to protect your clothing from the road while adding a sleek look to your chair.



SPECS

Transport Weight*	13.9 lbs.*	6.3 kg
Seat Width	14-22"	36-56 cm
Seat Depth	14-20"	36-51 cm
Front Seat Height	13.5-21"	34-53 cm
Rear Seat Height	11.5-20"	29-51 cm
Top Frame Tubing Size	1"	2.54 cm
Bottom Frame Tubing Size	1.125"	2.86 cm
Weight Capacity	300 lbs.	136.08 kg
Ansi Resna WC-19 Transit Approved	Yes	Yes

* 16"x16" Frame, 5"x1" Poly Casters, Height Adjustable T-Arm Receivers.
Footrest and Rear Wheels removed.



Stevens Point, WI 54481 | 715-254-0991 | www.kimobility.com

DCN0319.0