

New!

Get Out There!

GO-CHAIR

GET OUT THERE IN STYLE WITH THE ALL-NEW GO-CHAIR



NEW, dual under-seat compartments for convenient storage

Shown in: Ruby Red



The all-new Go-Chair is re-engineered from the ground up, offering a bold new style available in an array of contemporary colors. Enhanced performance and comfort, along with feather-touch disassembly, allows you to enjoy light-weight travel and independence on the go. With an increased weight capacity of 300 lbs. and a sleek, bold look, the new Go-Chair makes travel easy.



One-hand "Feather Touch" disassembly

Sapphire Blue Rose Quartz Citrine Yellow Amber Orange Pearl White



Variety of color options

SPECIFICATIONS

GO CHAIR


Drive Wheels	9" (solid)
Front Casters	6" (solid)
Rear Anti-tips	3" (solid)
Maximum Speed ¹	Up to 3.7 mph
Ground Clearance ²	1.6" at anti-tip bracket 1.9" at frame
Turning Radius ²	25.4" with foot platform
Length	33" with front rigging
Width ⁶	22.6"
Seat Type	Foldable, compact, lightweight
Seat to Ground Height	22"-23"
Batteries ^{3,5}	(2) 18 Ah, 12 V, deep-cycle
Range ^{1,4}	Up to 8.74 miles
Charger	Off-board, 2 amp
Electronics	40 amp, PG GC3 (non-programmable)
Weight Capacity	300 lbs.
Base Weight	67 lbs. (front & rear)
Standard Seat Weight	30 lbs.
Battery Pack Weight ⁵	31 lbs.
Suspension	Limited
Heaviest Piece	36 lbs. (rear)

- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition tire type and condition. These specifications can be subject to a variance of +10%, -5%.
² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
³ AGM or Gel-Cell type required. *See "Batteries and Charging."
⁴ Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
⁵ Battery weight may vary based on manufacturer.
⁶ Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA and ISO 7176 series standards. This product is not coded by DME PDAC or does not meet criteria: K0899. All specifications subject to change without notice.

COMPATIBLE LIFTS

- Outlander Full Platform
- Olympian
- Backpacker Series
- Commander

 To find a compatible lift for your vehicle/mobility product, use our new tool: pridemobility.com/lift-finder/

4-wheel

FEATURES

- Choice of 6 colors
- 18" x 17", 5-point swivel seat with latch release
- Increased 300 lb. weight capacity
- Faster speed - up to 3.7 mph
- Larger, 18 amp batteries
- Improved agility - longer wheel base; wider stance; compact turning radius 25.4"
- Dual swing-away, under seat storage bins
- Larger, height & width adjustable armrests
- Larger footplate
- One hand feather-touch disassembly
- 60" lap belt

OPTIONS

- Cup Holder
- Lap Belt- 50", 70"
- Safety Flag
- Oxygen Holder
- Dual Saddle Bag
- Weather Cover, Micro
- Additional Battery Pack
- Rear Mount Bracket
Cane/Crutch Holder
Single (3"),
Double (3")
- Armrest Mount Cane Clip

Pride
Mobility Products Co.
— LIVE YOUR BEST! —

GO GO
TRAVEL MOBILITY

Beamsville, Ontario · 888-570-1113

[PrideMobility](#) [PrideMobilityProductsCorp](#) [PrideMobilityProductsCanada](#)

pridemobility.com