

**NEW**



# Defender Reha CANADA

The Special Needs 2-in-1 CarSeat Solution  
Helping Families Enjoy Safe Travelling

Thomashilfen & **harmony** are two family owned companies passionate about providing safe and comfortable products. We are committed to developing and delivering products that help families enjoy safely travelling together.

## Let us **move** you!

“The Defender Reha allows children with special needs 360° of impact protection. The head, trunk and thighs supports, combined with our unique accessories, offer many positioning options. The seat allows for growth while remaining narrow enough for a second car seat as well as another child to fit in the back seat.”

Thomashilfen

Ultra plush seat pad made with durable premium fabrics and extra foam for maximum comfort is removable for machine washing



**Lateral trunk supports** for an upright, stable sitting position

9-position, 1-hand height adjustable headrest



**Head support** holds the child's head safely in place

Deep, supportive headrest, backrest and seat base for 360° impact protection



Chest clip and soft harness covers included



**Lateral side supports** are removable to adjust seat width for plenty of room for the child to grow with the seat



**SWIVEL BASE** option for easier transfer in **BOOSTER MODE**



Cup/snack holder swivels and can easily be installed on either armrest

The no-threaded 5-point harness provides additional support in Booster Mode for children weighing 18,1 to 50 kg / 40 to 110 lbs. Small crotch pad included.

Positioning-options **INCLUSIVE**



Smooth up-front harness adjuster

## Defender Reha CANADA

The Positioning 2-in-1 Special Needs CarSeat

## Accessories & Measurements



Swivel base with footrest adapter



UAS connector and stabilizing bow for swivel base



Footrest, short with footrest adapter and footrest pad



Footrest, long with footrest adapter



Seat wedge, below (+ 10°)



Table = playtray and additional impact shield



Abduction block



Crotch pad, large

Measurements & weights	Defender Reha CANADA
Seat depth	33 cm / 13.5"
Seat width*	21,5 - 30,5 cm / 8.5 - 12"
Seat to head**	60 - 78 cm / 23.5 - 30.5"
Shoulder height (child)**	25 - 49,5 cm / 10 - 19.5"
Lower leg length	(18 cm ***) 28 - 43 cm / (7" ***) 11 - 17"
Overall dimensions (w x h x d)	46,3 x 67,3 x 50,8 cm / 18.25 x 26.5 x 20"
Weight of seat	10 kg / 22 lb

### HARNES MODE

Weight of child	10 - 29,5 kg / 22 - 65 lb
Height of child	68,6 - 144,8 cm / 27 - 57"

### BOOSTER MODE

Weight of child	18,1 - 50 kg / 40 - 110 lb
Height of child	86,4 - 144,8 cm / 34 - 57"

- \* To reach max. seat width remove additional hip guide pads from additional seat pad
- \*\* To reach max. back height and shoulder height remove additional seat pad and foam pad under seat cushion
- \*\*\* To reach min. lower leg length you need to use footrest pad 834

Product & accessories	item code
<b>Defender Reha CANADA 2-in-1 CarSeat</b>	<b>888</b>
<b>Defender Reha CANADA BOOSTER Seat</b>	<b>887</b>
<i>(Child Weight 18,1 - 50 kg / 40 - 110 lbs)***</i>	
Swivel base with footrest adapter	820
UAS connector (LATCH) (for swivel base)	819
Stabilizing bow (for swivel base)	818
Footrest, short (depth 12 cm / 4.7")* (adjustable in angle)	806
Footrest pad (height 10 cm / 4") - (for 806)	834
Footrest, long (depth 18 cm / 7.1")* (adjustable in angle)	805
Footrest adapter + Base plate*	821
Seat wedge, below (+ 10°), black	810
Table, black	831
Crotch pad, large, black**	832
Abduction block, black	833

- \* The footrest adapter must be ordered to use footrest
- \*\* We recommend using larger crotch pad only in booster mode.
- \*\*\* Booster Mode Seat without UAS and top tether

Defender Reha **AND** all accessories  
exceed Canadian CMVSS 213 and  
US FMVSS 213 standards

