

INVACARE® MATRX®

# ELITE

BACK

---



---

The **Invacare® Matrx® Elite Back** is designed to provide precise orientation within the wheelchair for optimal postural support. The attractive, lightweight design appeals to individuals with active lifestyles.

The **Invacare® Matrx® Elite Full Back** with SuperSoft foam is designed for highly sensitive users or those who require accommodation of bony prominences. Improved immersion and pressure re-distribution for tilt/recline applications.

---

## Features and Functional Benefits

**Outer Cover** - Oversized cover/foam for extra comfort and skin protection. Moisture resistant and breathable. Optional privacy flap available.

**Hardware** - "Easy Set" hardware installs and adjusts in minutes – adjustments can be made while individual is seated in the wheelchair. Refer to pages 14/15 for additional hardware options available upon request.

**Back Shell** - Stylish, breathable and lightweight back shell. Heavy duty back shell includes extra layer of reinforcement. Shells are pre-drilled for headrest adapter attachment.

**Adjustable Lumbo-Sacral Support** - Includes removable positioning pad.

**Headrest Adapter** - Mounts directly onto Invacare Matrx Elite Back. Pre-drilled to accept headrest hardware. Available upon request.

**Optional Lateral Supports - Fixed** - 2.5" of vertical adjustment on back pan. Mounts inside back shell without interfering with depth/angle adjustment. Available in standard 5" H x 7.5" L and pediatric 3" H x 5.5" L sizes.

**Offset Fixed** - Provides an additional 3" of vertical position (higher or lower) relative to mounting slots. Available 5" H x 7.5" L.

**Swing-Away** - 2.5" of vertical adjustment on back pan. Mounts inside back shell without interfering w/depth/angle adjustment. Easy push button release, 4 pad sizes (see below).

**Offset Swing-Away** - Provides an additional 3" of vertical position (higher or lower). Mounts inside back shell without interfering with depth/angle adjustment. Easy push button release, 4 pad sizes (see below).

**Lateral pads For Swing away laterals - XS:** 3.25" H x 4.25" L, **S:** 3.5" H x 5.25" L, **M:** 3.75" H x 5.75" L, **L:** 4.25" H x 6.5" L



1. 3" Contour Depth
2. SuperSoft Foam Option
3. Fixed Lateral and Offset Fixed Lateral Pads
4. Swing Away
5. Offset Swing Away
6. Headrest Adapter Available



Model #	Height	Width	Depth Adjustment	Angle Adjustment	Product Weight	Weight Capacity	Limited Warranty
PBE	10", 12", 14"** 16", 18", 20"***	12" - 20"	4" (2" in front or behind back canes)	40° range +/- 20°	1.7 lb. and up	300 lb.	2 Years
PBE-FB**	20"	16", 17", 18", 19", 20"	4" (2" in front or behind back canes)	40° range +/- 20°	1.8 lb. and up	300 lb.	2 Years
+PBE- HD	16", 20"	21"-24" (25"-30" Available Upon Specials Request)	4" (2" in front or behind back canes)	40° range +/- 20°	6.1 lb. and up	500 lb.	2 Years

\* 14" high now available in 12"-18" widths

\*\* Fullback includes SuperSoft foam - laminate of 1" firm, 1" soft HR foam for improved immersion and comfort