




# GoTo

**PRESSURE  
RELIEF CUSHION**




# KEYPOINTS

Combination of memory and support foams provide moderate pressure relief and improved pressure distribution to GoTo Seat users.






GoTo without pressure relief cushion



- 90.0
- 84.3
- 78.7
- 73.0
- 67.3
- 61.7
- 56.0
- 50.3

GoTo with pressure relief cushion



- 44.7
- 39.0
- 33.3
- 27.7
- 22.0
- 16.3
- 10.7
- 5.0
- 0.0

Scale and set-up remained constant.

# KEYPOINTS

---

The Cushion and Roll also provide postural support and pelvic stability – placing the roll underneath the contoured cushion will aid comfort and stability and reduce stretch on gastrocnemius (a calf muscle) for those in corrective footwear or AFOs.



# OTHER WAYS TO USE THE CUSHION AND ROLL

---



**Long sitting in the floorsitter**

**On a dining  
room chair**



# GENERAL INFORMATION

---

## SIZE 1

167-101



## SIZE 2

167-102



---

Can be purchased with GoTo or as an accessory.  
Accessory – **£59** / Purchased with GoTo – **£55**

# FEATURES & BENEFITS

Useful at school,  
when out and about,  
or watching TV and  
relaxing at home

Contoured cushion provides additional pressure relieving comfort as well as pelvic stability and postural support, whilst roll can be used in a variety of different ways to enhance secure and comfortable positioning

Cushion and roll can be used together or separately – attached by Velcro

Roll also helps to prevent rubbing or friction (legs against hard chair surfaces)

Non-slip base

Easily attached to GoTo Seat by Velcro strap

