

INVACARE® MATRX®

VI CUSHION



The **Invacare® Matrx® Vi Cushion** is designed to provide skin protection and positioning for individuals at moderate to high risk of skin breakdown. The unique three part foam construction targets skin protection where it is needed most, and features excellent stability, support and comfort.

Features and Functional Benefits

Outer cover - Moisture resistant and breathable. Includes slip resistant base, hook and loop fasteners, lifting strap. Available in black and Kid*ab*ra print.

Inner cover - Moisture resistant. Easily wipes clean. Zippered cover allows access for cushion modifications.

Top layer - Soft (ILD) High Resilient (HR) foam for optimal comfort and support. Pre-ischial contour provides stability to pelvis. Abductor/adductor contour promotes midline leg positioning. Made with Ultra-Fresh™ for antibacterial and odor protection.

Visco foam insert - Molded visco-elastic foam ischial relief insert. Anatomically targeted to ischial/sacral area. Provides immersion to sensitive bony prominences.

Bottom layer - Firm (ILD) High Resilient (HR) foam for support and durability. Trochanteric shelf provides stability to hips. Waterfall front edge. Rear radius for sacral support.

Drop Base - Includes built-in drop base, drop hooks, DuraClip fasteners. Allows adjustment of seat angle.

Pelvic Obliquity Kit - Easily installed on all Invacare Matrix Vi/Kid*ab*ra-Vi cushions/bases. Accommodates 0.5" or 1" obliquities.

Custom Modifications - Leg length discrepancies, width/depth reductions. Foam, inner and outer covers modified to your specifications.

Cushion Ridigizer - Lightweight, slim profile solid seat insert. Contoured to eliminate "hammock effect" of sling upholstery.



1. Drop Base
2. Kid*ab*ra Fabric Standard on Pediatric Sizes
3. Custom Modifications
4. Cushion Ridigizer
5. Pelvic Obliquity Kit
6. Rear Cushion Radius Diagram



Model #	Width*	Depth**	Height at Adductor/Abductor	Height at Leg Trough	Product Weight	Weight Capacity	Limited Warranty
KA-VI	12" - 14"	12" - 14"	2.87"/2.75"	2"	1.5 lb. and up	150 lb.	2 Years
MA(WWDD)-VI (Standard)	14", 16", 18", 20", 22"	15" - 20"	3.75"/3.25"	2.5"	3 lb. and up	300 lb.	2 Years
MA(WWDD)-VI (Heavy Duty)	18" - 30"	18" - 26"	3.75"/3.25"	2.5"	4 lb. and up	600 lb.	2 Years

* Actual width is 1/2" narrower than stated width

** Actual depth is 1 1/2" longer than stated depth