

INVACARE® MATRX®

**PB**  
**BACK**

---



---

The **Invacare® Matrx® PB Back** is contoured to provide midline support and positioning. The extensive range of sizes provides exceptional support and comfort for a variety of individuals.

---

## Features and Functional Benefits

**Foam** - Molded foam made with Ultra-Fresh™ for antibacterial and odor protection wraps around shell for maximum comfort and protection. Foam attaches to back shell via hook and loop strips for easy replacement.

**Outer Cover** - Moisture resistant and breathable. Available in black or Kid\*ab\*ra print.

**Hardware** - "Easy Set" hardware is easily installed and adjusted in minutes – adjustments can be made while individual is seated in the wheelchair. Refer to pages 14-15 for additional hardware options (Available upon request).

**Back Shell** - Contoured shell provides exceptional midline support and positioning. Durable, lightweight ABS back shell

includes built-in plate for headrest attachment. Heavy-Duty sizes (20" wide and up) include an aluminum back shell. 500 lb. weight capacity

**Visco Foam Upgrade** - Laminate of 0.5" high resilient foam and 0.625" Visco Foam. Available on Invacare Matrx PB and Elite Backs. Molded foam is not available on Invacare Matrx PB Back with this option.

**Removable Lumbo-Sacral Support** - Includes lumbo-sacral support that is easily adjustable to meet individual's needs.

**Foam-in-Place Kit** - Individually molded insert produced in minutes with user seated in back. Optional high stretch cover available.

**Privacy Flap** - Optional privacy flap is available.



1. 5" Contour Depth
2. Kid\*ab\*ra Fabric Optional
3. Privacy Flap
4. Foam-in-Place Kit\*
5. Visco Foam Upgrade
6. Adjustable Sacral/Lumbar Support



Model #	Height	Width	Depth Adjustment	Angle Adjustment	Product Weight	Weight Capacity	Limited Warranty
MPB (Standard)	12", 16", 20"	14", 16", 17", 18", 20"	4*** (2" in front or behind back canes)	40° range**	2.8 lb. and up	300 lb.	2 Years
HDB (Heavy-Duty)	16", 20"	20", 22", 24"	4*** (2" in front or behind back canes)	40° range**	6 lb. and up	500 lb.	2 Years
HDBxW16 (Extra Wide) HDBxW20 (Extra Wide)	16", 20"	21", 23", 25"-30"	4*** (2" in front or behind back canes)	40° range**	6 lb. and up	500 lb.	2 Years

\*High stretch cover included when backs are ordered with Foam-in-Place option

\*\* Limitation of depth/angle adjustment if used on chairs narrower than model width