



Shown in Merino  
Fawn Fabric



Shown in Badlands  
Walnut Fabric

# *Viva* LIFT!™

## POWER RECLINERS

By Pride® Mobility Products Corp.

### ***Breakthrough Technology!***

**Power Lumbar** for personalized ergonomic comfort

**Full-Width Power Headrest** to lift head,  
neck and shoulders

**Lithium Battery Back-Up** is an industry first!  
Abundant power lasts for days and replaces  
ordinary one-time-use backups

**Fashion Forward Styles** set these chairs  
apart from all others!

Elevate Your **Lifestyle!**

Beamsville, Ontario | 888-570-1113  
pridemobility.com | prideliftchairs.com

**Pride**®  
Mobility Products Co.  
— LIVE YOUR BEST —

### Elegance



STANDARD FEATURES:

- POWER HEADREST
- POWER LUMBAR
- USB REMOTE
- LITHIUM BATTERY BACKUP
- ALL CLOSE AND LIFT
- SPRING LOADED OTTOMAN
- INFINITE LAY FLAT

### Urbana



STANDARD FEATURES:

- POWER HEADREST
- POWER LUMBAR
- USB REMOTE
- LITHIUM BATTERY BACKUP
- ALL CLOSE AND LIFT
- SPRING LOADED OTTOMAN
- INFINITE LAY FLAT

### Perfecta



STANDARD FEATURES:

- POWER HEADREST
- POWER LUMBAR
- USB REMOTE
- LITHIUM BATTERY BACKUP
- ALL CLOSE AND LIFT
- SPRING LOADED OTTOMAN
- INFINITE LAY FLAT

### Metro



STANDARD FEATURES:

- POWER HEADREST
- POWER LUMBAR
- USB REMOTE
- LITHIUM BATTERY BACKUP
- ALL CLOSE AND LIFT
- SPRING LOADED OTTOMAN
- INFINITE LAY FLAT

Model	PLR975M	PLR965M	PLR945M	PLR925M
Position	Relaxed Reclined Infinite Lay Flat Upward Lift 5'4" - 6'	Relaxed Reclined Infinite Lay Flat Upward Lift 5'4" - 6'	Relaxed Reclined Infinite Lay Flat Upward Lift 5'4" - 6'	Relaxed Reclined Infinite Lay Flat Upward Lift 5'4" - 6'
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	20"	20"	19"	19"
Seat Depth	20"	20"	20"	21"
Seat Width	22"	21"	22"	21"
Top of Back to Seat	27"	25"	25"	26.5"
Back Style	Blown fiber, Pillow back	Blown fiber, Split back	Blown fiber, Split button back	Blown fiber, Split back
Fully Padded Chaise	Yes	Yes	Yes	Yes
Footrest Extension	Standard	Standard	Standard	Standard
Distance From Wall	14"	14"	15"	14"
Height	43"	42"	43"	41"
Overall Width	35.5"	33"	31"	33"
Chair Weight	152 lbs.	145 lbs.	142 lbs.	142 lbs.
Warranty	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

STANDARD FABRIC SELECTION GUIDE

BADLANDS:

Mushroom Walnut

STONEWASH:

Granite Gun Metal

MERINO:

Fawn Steel

SAVILLE:

Brown Grey

POSITIONS:

