



# POWER PLUS MOBILITY

**BUILT FOR YOU**  
TECHNOLOGY



**STP**

*For more information or pricing please contact  
toll free 1.855.614.0333 or visit our website @ [powerplusmobility.com](http://powerplusmobility.com)*

## FEATURES

- Light Weight Frame
- Depth Adjustable Frame
- ABS Seat Pan
- Dual Gas Struts for Easy Tilt (up to 45 degrees)
- New Recessed Lower Frame For Foot Propulsion
- Dual Gas Cylinders for Smooth & Easy Tilting
- Relax Seating System
- New User Friendly Back Rest Hardware
- New Gas Cylinder Release Mechanism
- Silver Grey or Black Frame Colour
- T-Style Arms
- Adjustable Stroller Handle
- Swing In & Out Footrest
- Adjustable Caster Journals
- Recessed Axle Plates



**ADP CODE: WM0001040**

## TILT

- Infinite Tilt Adjustments
- Dual Gas Struts for Easy Tilt (up to 45 degrees)
- Two-Way Locking Gas Strut
- Single Release, Dual Cylinders
- Two Position Cylinder Mount

## UPCHARGE FEATURES:

- Seat Width Over 18"
- Seat Depth Of 14", 15" & Over 18"
- Custom Back Height
- Laterals
- Custom Headrest
- Upgrade Seating (Relax Seating System)
- Elevating Legrest
- One Piece Elevating Legrest System
- Angle Adjustable Footplate
- One Piece Footboard
- Extra Large Footplate
- Contracture Sling
- Non Standard Arm Height
- Recline Back
- Non-Standard Frame Colours

## STANDARD SPECS.

Seat Width:	14" - 18"
Seat Depth (Adjustable):	16" - 18"
Seat Height (Adjustable):	12.25" - 18.25"
Seat Angle(Infinite):	0° - 45°
Back Angle (Adjustable):	0°, 5°, 10°, 15°
Adjustable Stroller Handle:	0° - 180°
Weight Capacity:	250 lbs
Weight of Chair:	60 lbs (approx.)

## SEE ORDER FORM FOR CUSTOMIZATION

\*SEE ORDER FORM

Warranty:

5 Yrs Frame, 2 Yrs Components

**"Our Products Will Move You"**

**POWER PLUS MOBILITY**

Tel: 905.614.0333 • Fax: 905.614.0435  
208 Wilkinson Road, Brampton L6T4M4



website: [www.powerplusbility.com](http://www.powerplusbility.com)