



X'CAPE™

The Zippie X'CAPE weighs as little as 19 pounds and is the first ever wheelchair with XLOCK™ Technology for the performance of a rigid frame and convenience of a folding frame.



Lily explores her world in the Zippie X'CAPE

FEATURES AND BENEFITS

Patented XLOCK™ Technology

XLOCK lockout crossbrace and fewer moving parts reduce weight and frame play for greater efficiency and less fatigue.

Built-in Growth

2" of width growth in the crossbrace and 3" of depth growth in the seat rail keep your child positioned properly as they grow.

WC-19 Transit Approved

WC-19 transit approved for use as an occupant's seat in a motor vehicle when used with an approved Wheelchair Tiedown and Occupant Restraint System in accordance with SAE (J2249).

Z-finity™ System

Lower leg length setting can be adjusted between -.5" and 21" from the seat and footrest angle can be adjusted infinitely to achieve the best positioning as your child grows.



Swing-in/Swing-out & Fixed Hanger Options

Swing-in/Swing-out hanger option aids in transfers and pulling up close to cabinets and counters. Fixed front frame option offers the lightest weight configuration and an active style.

VersaRail™ Seat Rail

VersaRail provides your child access to all of the components and accessories needed – custom seating, armrests, wheel locks, side guards, etc. – even on the shortest frames.



One-step Push-to-Lock
One-step push-to-lock tab locks out the frame.



One-step Pull-to-Fold
One-step pull-to-fold strap collapses the frame effortlessly.



Folding Push Handle
Tucks away or raises to your height.



Color Options

27 frame colors and 24 crossbrace colors can create over 600 combinations!

Refer to the Quickie®/Zippie® Color Brochure or the Sunrise Medical website to view color choices.

HCPCS Code: E1238

Chair Weight: 19 lbs

8x8" fixed front frame, 20" lite spoke wheels. Excludes armrests, push handles & anti-tips. Actual weight may vary according to chair configuration & accessories.

Overall Width
14.8" to 26.5"

Folded Width
Starting at 5" (without wheels, hangers, and footplates)

Weight Capacity
165 lbs (w/o seating)

Seat Widths
8" to 16" (grows to 18")

Seat Depths
8" to 18" (grows to 20")

Front Seat-to-Floor Heights
13" to 19"

Rear Seat-to-Floor Heights
11" to 18.5"

Footrest/Frame Angle
80° (Zfinity™ may be adjusted to achieve a wide range of footrest angles)

Transit Option: Yes
Conforms with ANSI/RESNA WC-19 standards. Refer to Quickie® & Zippie® WC-19 Transit Securement Supplement (MK-100072) for additional information. Supplement is available for download at www.sunrisemedical.com.



800.333.4000 www.SunriseMedical.com

©01.2015 Sunrise Medical (US) LLC
MK-100118 Rev.A

Aztec Gold	Black	Black Cherry	Black Opal
Blue Opal	Candy Blue	Candy Purple	Candy Red
Electric Blue	Evergreen	Glow	Green Apple
Hot Sparkle Pink	Mauve Pink	Rootbeer	Sparkle Silver
Sunrise Orange	Yellow	Titanium Color Paint	Matte Black
Matte Black Cherry	Matte Purple	Matte Electric Blue	Matte Evergreen
Mossy Oak Camo*	Stars and Stripes*	Zebra*	*Kolorfusion not avail. on the crossbrace

Click here for printer-friendly version