

## Product description

Our T350 model is an innovative and **special needs trike**. Just like its little brother the T-250, it can be easily adapted by simply removing, adding or adjusting features, without any tools! The T-350 has the longest life cycle on the market, as it can adapt to a child from 10 years old up into adulthood. With a



total width of 30", the T-350 can easily fit through any standard doorframe. It can thus be used at school, in a therapy

centre or at home. In addition, its simple disassembly mechanisms allow you to store it in your car. Designed mainly to help kids with special needs to strengthen their muscles and improve their coordination and balancing skills, the T-350 is also a recreational tool which gives them the chance they deserve to ride, have fun and interact socially with friends and family. You can pull the trike with an optional tow bar, if required, thanks to a simple mechanism allowing you to switch between fixed or freewheel gearing systems. The T-350 can be set at 6 or 18 gears. There are more than 20 adjustable and removable options to suit any child's specific physical needs. It is available in four colour choices.

## Features and options

- Adjustable mesh backrest
- Adjustable and removable headrest
- Removable chest harness
- Speed reducer
- Drilled pedal crank
- Parking brakes



- Abduction bar
- Hip and thoracic seatbelt and additional support pads
- Mustang handlebars
- Direction spring
- Removable abduction pommel
- Adjustable pedals
- Vertical and horizontal seating adjustment
- Fixed or freewheel gearing system
- And more!



## Specifications

Bike weight	77 lbs
Total length	70"
Frame height	42"
Wheel size	20"
Width	30"

## Sizing details

User age	10+ years old
Leg length (measurement from greater trochanter to floor)	26" to 45"
Maximum load	250 lbs

