

QUICKIE

CANADIAN-MADE POWER WHEELCHAIRS



The Quickie Power X series is a family of “Made in Canada” power wheelchairs providing reliable mobility to thousands of Canadians since 2006.

NEW GENERATION

X series Power Wheelchairs



The new generation offers many options including a new shroud to give your power wheelchair a new look and style. When performance meets beauty the perfect power chair is waiting for you.

What is unique to the X series?

Smooth responsive ride

No adjustments needed to achieve stability and traction.

Top of the line comfort

While riding over uneven surfaces, the hard movements and vibrations are controlled and not absorbed by the seat or user.

Worry free travel over obstacles

The unique suspension will distribute the weight over the entire footprint.

Intellidrive

A unique steer correcting feature that improves maneuverability in indoor or outdoor environments allowing you to drive with ease - naturally.

Xperience 2™

MID-WHEEL



Recline



Tilt



New Shroud

The Quickie® Xperience 2™ is a mid-wheel drive power wheelchair offering excellent stability for better comfort. A short turning radius gives easy access to small spaces. The independent suspension provides users with a smooth ride by dampening the effect of day-to-day obstacles. The Xperience 2™ is one of the best options for users spending time inside and out. It is available in Standard, Heavy Duty or Pediatric.

SPECIFICATIONS

Weight Capacity	300 lbs. Standard, 450 lbs. Heavy Duty option
Front Caster	6" Caster, 8" Caster pneumatic, with or without insert
Drive Wheel	12.5" and 14" Cast wheel brackets, pneumatic with or without insert
Transit	4-Point wheelchair transport bracket included
Electronics	1. VR2 90 Amp standard non expandable 2. R-Net EL 90 Amp expandable, LED joystick 3. R-Net 120 Amp expandable, LCD screen joystick
Seat to floor height	16" Note: Can be as low as 15" with 12.5" drive wheels
Max Speed	6.5 MPH
Overall Base Width	25.5" (Can be as narrow as 23.5" with 12.5" drive wheels and 22NF battery box)
Overall Base Length	43.5"
Seating	Fully adjustable seat: Rehab Seat
Warranty	5 years for frame, please see owners manual for all warranty information
Funding Info ADP#s	PEDIATRIC: WPK3QU080 ADULT: WPA3QU280 Under the AADL program: W699

Colours available

Candy Red	Black Opal Metallic
Midnight Purple	Razzleberry
Matte Black	Black
Pearl Pink	Toxic Green
Silver Metallic	Candy Blue
Yellow	Burnt Orange Metallic
Gun Metal	These colours are for a guideline only. Please see dealers for true colour match.
Sparkle Orange	

Xplore 2™ HYBRID

REAR WHEEL



New Shroud



New Wheel option



New Back Cover

The Quickie® Xplore 2™ is the most hybrid product available. This wheelchair offers the best of both worlds; the traction of a rear wheel drive and the short turning radius of a mid-wheel drive wheelchair. The fast speed of the Xplore 2™ makes it the ideal choice for someone spending a lot of time outdoors and who wants to go the distance. Safe and comfortable, the 6-wheel independent suspension is offered in Standard and Heavy Duty. The wheelchair come standard with transport brackets. Seat options: tilt, manual and power recline, and elevating seat.

SPECIFICATIONS

Weight Capacity	300 lbs. Standard, 400 lbs. Heavy Duty option
Front Caster	8" Caster, 9" Caster pneumatic, with or without insert
Drive Wheel	12.5" and 14" Cast wheel pneumatic, with or without insert
Transit	4-Point wheelchair transport brackets included
Electronics	1. VR2 90 Amp standard non expandable 2. R-Net EL 90 Amp fully expandable controller, LED joystick 3. R-Net 120 Amp fully expandable controller, LCD Screen Joystick
Seat to floor height	16"
Max Speed	8.5 MPH
Overall Base Width	25.5" (As narrow as 23.5" with 12.5" drive wheels and 22NF battery box)
Overall Base Length	43.5"
Seating	Fully adjustable seat
Warranty	5 years for frame, please see owners manual for all warranty information
Funding Info ADP#s	PEDIATRIC : WPK3QU070 ADULT: WPA3QU270

Colours available

Candy Red	Black Opal Metallic
Midnight Purple	Razzleberry
Matte Black	Black
Pearl Pink	Toxic Green
Silver Metallic	Candy Blue
Yellow	Burnt Orange Metallic
Gun Metal	These colours are for guidelines only. Please see dealers for true colour match.
Sparkle Orange	

XCEL 2™

BARIATRIC



Heavy Duty Tilt



The mid-wheel Heavy Duty Quickie® XCEL 2™ features the Advance Geometric Design (A.G.D) which helps enhance the stability, mobility and freedom for the end user. Engineered without compromise, the XCEL 2™ offers a weight capacity of 550 lbs., a bariatric CG power tilt bariatric power recline option, and heavy duty elevating legrests or center mount for the best comfort.

SPECIFICATIONS

Weight Capacity	550 lbs. Standard
Front Caster	8" Caster
Drive Wheel	14" Cast wheel pneumatic, with or without insert
Transit	4-Point wheelchair transport bracket included
Electronics	VR2 90 Amp drive controller R-Net EL 90 Amp expandable controller - LED hand control R-Net PM 120 amp fully expandable controller
Seat to floor height	16.5" - with standard seating 17 1/4" - with tilt
Max Speed	5.4 MPH
Overall Base Width	27 3/8"
Overall Base Length	43.5"
Seating	HD fully adjustable seat, bariatric Eclipse
Warranty	5 years for frame, please see owners manual for all warranty information

Colours available

Candy Red	Black Opal Metallic
Midnight Purple	Razzleberry
Matte Black	Black
Pearl Pink	Toxic Green
Silver Metallic	Candy Blue
Yellow	Burnt Orange Metallic
Gun Metal	
Sparkle Orange	

These colours are for a guideline only. Please see dealers for true colour match.

EVERYDAY OPTIONS



Infinite adjustment
center mount -5° to 70°



Light Package



New cantilever armrest
opens with 5° angle



New Back Cover



Phone & Tablet Support



Independent height
adjustable footplates



Coffee Cup Holder



Colour Headrest Mount



New Drive Wheel

Please see order forms for a complete list of options



zippie® Xseries

zippie® Xperience 2™



zippie® Xplore 2™ HYBRID



The Zippie Xperience 2 is the wheelchair that is versatile whether you are inside or at the park! The mid-wheel pediatric wheelchair offers a fully adjustable rehab seat, stability, maneuverability and provides independence to let kids be kids!

Weight Capacity	300 lbs. Standard
Front Caster	8" Caster, 9" Caster pneumatic, with or without insert
Drive Wheel	12.5" and 14" Cast wheel pneumatic, with or without insert
Transit	4-Point wheelchair transport bracket included
Electronics	1. VR2 90 Amp standard non expandable 2. R-Net EL 90 Amp fully expandable controller, LED Joystick 3. R-Net 120 Amp fully expandable controller, LCD Screen Joystick
Seat to floor height	15 3/4"
Max Speed	6.5 MPH
Overall Base Width	25.5"
Overall Base Length	43.5"
Seating	Fully adjustable seat
Warranty	5 years for frame, please see owners manual for all warranty information
ADP#	PEDIATRIC : WPK3QU080

For kids who want to go fast, and do not want to be stopped by any obstacles out there, the Zippie Xplore Hybrid 2 is the right choice. Its high speed, short turning radius and multiple configurations make it the best option for young kids and teenagers who are busy with everyday activities.

Weight Capacity	300 lbs. Standard
Front Caster	8" Caster, 9" Caster pneumatic, with or without insert
Drive Wheel	12.5" and 14" Cast wheel pneumatic, with or without insert
Transit	4-Point wheelchair transport bracket included
Electronics	1. VR2 90 Amp standard non expandable 2. R-Net EL 90 Amp fully expandable controller, LED Joystick 3. R-Net 120 Amp fully expandable controller, LCD Screen Joystick
Seat to floor height	16"
Max Speed	8.5 MPH
Overall Base Width	25.5"
Overall Base Length	43.5"
Seating	Fully adjustable seat
Warranty	5 years for frame, please see owners manual for all warranty information
ADP#	PEDIATRIC : WPK3QU070





Complement your wheelchair with BACKS & CUSHIONS



The JAY® J3™ cushion provides outstanding stability and optimal pressure distribution for the client at high risk of skin breakdown, and offers the choice of JAY Flow™ Fluid or ROHO DRY FLOATATION® air inserts.



The JAY® Union™ is a versatile, comfortable, skin protection and positioning cushion composed of a dynamic fluid and foam layering system, moisture resistant inner cover, and an anti-microbial outer cover.



The JAY® Fusion™ is an adjustable cushion designed for skin protection that features JAY Flow™ Fluid or ROHO DRY FLOATATION® air insert options in a package that's comfortable, stable and supportive.



The JAY® Ion™ is a comfortable skin protection and positioning cushion composed of a multi-layered foam base, moisture resistant inner cover, and an anti-microbial outer cover.



The JAY® J3™ backrest product line offers a multitude of width, height, and contour depth back shells to fit the back to the client. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.




Comfort and Support all in One! The JAY Union Back is the complementary addition to any manual, positioning on power wheelchair.



American Access

ACCESSIBILITY RAMPS

Our wide range of accessibility ramps make mobility even easier. They are easy, secure and flexible to use, along with a great warranty as well as colour options to accentuate any residence, and much more.



- Entrada Ramp Series
- Aluminum Threshold Ramps
- Plastic Threshold Ramps
- Hybrid Ramp Series
- XM Ramp Series
- Folding Sidekick Ramps



Customer Service: 1.905.660.2459 / 1.800.263.3390
or visit www.SunriseMedical.ca

OXPOWER SERIES