

Advantage Rail<sup>™</sup> / Portable Options A floor mounted pivoting and locking horizontal rail. Excellent where ceiling restrictions prevent a floor to ceiling pole. Perfect beside a toilet, bath, bed, or chair.

# Features and Benefits

#### Simple.

A gentle lift unlocks the rail to start it moving. Keep the rail up to move freely. Let the rail down to lock. Provides 8 locking positions. Height adjustable.

### Safe, Effective Support.

Push or pull on the rail to assist with sitting, standing and moving. To help prevent a forward fall, the rail will lock if released.

### Reduce Caregiver Strain.

Maximize the use of available strength while a caregiver supervises. Promotes safe and independent movement.

### Portable Model (sold separately).

Rail quickly secures and releases from a low profile, floor mounted plate. Install additional floor plates to use rail in multiple locations.





## **Technical Specifications:**

Advantage Rail™	17" / 43 cm length x 1.25" / 3.175 cm diameter grip. Locks every 45° (8 positions). Height range 30-38" / 75-95 cm cm from floor.
Grip Length	16" / 41 cm. Washable, latex free grip.
Advantage Rail™ Mast	28" / 71 cm height x 2.25" / 5.7 cm diameter.
Floor Plate	7.5" / 19 cm diameter floor plate (screw pattern 6 x 3/8" holes).
Weight Capacity	300 lbs / 147 kg. Bariatric weight capacity 450 lbs / 205 kg.

