

The Power EZ™ Pool Lift

Item# F-03PEZ

The Power EZ™ Lift is a great economical choice for individuals needing to gain access to their pool or spa.

The Power EZ™ Lift provides a fully automatic actuator to raise and lower with 360° manual rotation. The lift can accommodate spa walls up to 26 inches (66cm) high and can be set back up to 20 inches (50cm) from the pool edge.

The Power EZ™ Lift now comes with a rotation bar to allow attendant greater ease in rotating the unit.

Features:

- 400 Pound (181 kg) weight capacity
- Rechargeable battery operation with charger
- Clears walls up to 26" (66 cm) high
- Durable mesh sling- standard size: 38"x20"
- Adjustable stainless steel chains
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty

Options:

- Hard seat option Item #F-0320HSO
- Youth sling Item #F-34EZCS /34"x18" (86x46cm)
- Extra long sling Item #F-037EZLS / 50"x28" (127x71cm)
- Standard sling with rigid insert & seat belt -Item #F-036RSS
- Additional standard in ground anchor -Item #F-033EZA



Sling with rigid seat and belt



Hard Seat Option







- Wood deck anchor Item #F-30600
- Surface mount anchor Item #F-035EZSM
- Upgrade Pack Item #F-03PEZAP Includes battery, rigid sling with seat belt, sling bars and cover
- Cover Item #F-440EZC
- Hard seat cover Item #F-440EZHSC

PROUDLY SOLD BY: