

The Pro Spa 40 Lift™ Item# F-006PSLB

The Pro Spa 40 Lift[™] is a fully powered, ADA compliant* lift designed especially for use with above ground spas and pools with wall heights of 20" up to 40".

The lift provides easy access at the push of a button and is constructed of low carbon stainless steel and protected with a durable white epoxy powder coated finish.

Field measurements are required.



Clears walls up to 40" high!



- 400 Pound (181 kg) weight capacity
- ADA compliant for walls up to 37" high
- Clears walls up to 40" high, 14" thick
- Rechargeable battery operation with charger
- Submergible hand held remote
- Pull out leg rest
- Dual flip up arms
- Limited lifetime structural warranty

Options:

- Upgrade Pack Item #F-005SLBAP (must be purchased with lift) Includes extra battery, headrest, chest strap & cover
- Headrest Item #F-422PLH (Can't retrofit)
- Chest strap Item #F-423CS
- Wood deck anchor Item #F-106SLWDP
- Protective cover Item #F-445SLC





