



BED COMPARISON

Manufacturer	Model	Bed Width	Bed Length	Height	Power	Weight Capacity	Brakes	Surface	Reverse Trendelenburg
Carrol	Arro	35.5"	80"	6.75" - 33"	Head/Foot Height	500 lbs	Locking Casters	Pan	Yes
Carrol	Home Care	35"	80"	13.5" - 23.5"	Head/Foot Height	450 lbs	Floor Jack	Pan	No
Invacare	Echo	35"	76" or 80"	13.5" - 23.5"	Head/Foot Height	450 lbs	Floor Jack	Pan	No
Invacare	Full Electric	36"	80"	15" - 23"	Head/Foot Height	450 lbs	Locking Casters	Spring	No
Invacare	Full Electric Low	36"	80"	9.5" - 20"	Head/Foot Height	450 lbs	Locking Casters	Spring	No
Invacare	Manual	36"	80"	15" - 23"	All Manual	450 lbs	Locking Casters	Spring	No
Invacare	Semi-Electric	36"	80"	15" - 23"	Head/Foot Height	450 lbs	Locking Casters	Spring	No
Liftcare	Protean 3	35.5"	79" - 86"	3.5" - 31"	Head/Foot Height	550 lbs	Locking Casters	Pan	Yes
Rotec	Multiposition	35", 39" 54" or 60"	75" or 80"	13" - 23"	Head/Foot Height	250 lbs	Locking Directional Casters	Pan	Yes
Rotec	Multiposition 500	35", 39" 54" or 60"	75" or 80"	13" - 23"	Head/Foot Height	500 lbs	Locking Directional Casters	Pan	Yes
Rotec	Multitech	35", 39" or 48"	75" or 80"	12" - 30"	Head/Foot Height	500 lbs	Locking Casters	Frame	No
Rotec	Multitech 8	35"	78"	8.5" - 29.5"	Head/Foot Height	580 lbs	Central Brake	Frame	No