

Jazzy®

ELITE ES PORTABLE

The Jazzy® Elite ES Portable delivers a reliable blend of power, performance and portability. In-line, front-wheel drive technology gives the Elite ES Portable excellent stability and maneuverability while simple disassembly makes it easy to take along with you.



Jazzy® Elite ES
Portable shown
disassembled

Pride[®]
Mobility Products Co.

Specifications

Weight capacity	300 lbs.
Maximum speed	Up to 4 mph
Ground clearance	1.5"
Turning radius	25.5"
Overall length	34.75"
Overall width	22.75"
Drive wheels	9" solid
Front wheels	3" solid anti-tips
Rear wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC 3 controller
Specialty controls	N/A
Battery charger	Off-board, 3 amp
Per-charge range	Up to 11 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight	24.5 lbs. each
Base weight	Front: 47.6 lbs. Rear: 21 lbs.
Standard Seat weight	43 lbs.
Max. seat dimensions	20" x 20"
Warranty	5-year limited warranty on frame; 13-month limited warranty on electronics; 2-year limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

Features

- Easy disassembly for convenient transport
- In-line, front-wheel drive technology
- 40 amp, PG GC 3 controller
- Simple main frame design for easy serviceability

Options

- High-back seat with solid seat pan
- Depth-adjustable, limited recline high-back
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

Recommended Lifts

- Backpacker® Plus
- Outlander
- Backpacker® MV

Color Choices



Jazzy
Red



Jazzy
Blue