

HIGH-PERFORMANCE OPERATION AND ALL NEW FEATURES AT A GREATER VALUE THAN EVER BEFORE

VICTORY® TWIN

The sleek, sporty **Victory**® **Twin** scooter delivers high-performance operation, all new features and Feather-touch disassembly. With a full complement of advanced standard features like exclusive low-profile tires, dual front wheels for added stability and wraparound delta tiller, the Victory® Twin is a greater value than ever before.



Mobility Products Co.

Beamsville, Ontario · 888-570-1113 · www.pridemobility.com

Victory® Twin	
Model number:	SC712
Weight capacity:	350 lbs.
Maximum speed1:	Up to 5.25 mph
Ground clearance ² :	2.25" at motor
Turning radius:	42"
Overall length ² :	46"
Overall width2:	22.25" (at wheels)
Range per charge ^{1,4} :	Up to 15.3 miles
Tire size:	Front: Dual 3.5 x 9 solid Rear: 3.5 x 10 solid
Total weight without batteries:	149 lbs.
Heaviest piece when disassembled:	Front: 53.5 lbs.
Standard seating:	Type: Foldable molded plastic Weight: 27.5 lbs. Material: Black Vinyl Dimensions: 18" width 17" depth
Seat-to-deck range:	17.5" - 19"
Drive system:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual braking system:	Regenerative and electromechanical
Battery requirements ³ :	Two 12-volt, deep cycle Size: 40 AH or U-1 Weight: 40AH: 32.5 lbs. each U-1:24.5 lbs. each
Battery charger:	3-amp, off-board (U-1) 5-amp, off-board (40 AH)

Colors:



¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition.

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³AGM or Gel-Cell type required.

⁴Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.

 $NOTE: All\ specifications\ subject\ to\ change\ without\ notice.$

Features

- Dual front wheels for added stability
- Full lighting package, including high visibility automatic rear brake light system and turn signals
- Exclusive, one-hand Feather-touch disassembly provides the ultimate quick transport system
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Pride's exclusive stylish, light weight, non-marking, low profile wheels
- · Backlit battery gauge
- Long lasting, high-intensity LED headlight offers optimal pathway illumination
- Lightweight seat with viscoelastic foam inserts for increased comfort
- Seat includes pinchless hinges & 8-position seat rotation and sliders
- Easy-grip tiller adjustment knob
- Easy access, tiller-mounted charger port for standard dual voltage offboard charger
- · 4-Pole high-speed motors standard

Options

- High-back seat
- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- · Rear basket
- Safety flag
- Walker holder
- · Weather cover
- Saddlebag

Pride Lifts & Ramps





The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.

