

Dependa-Bar<sup>™</sup>
Dependa-Bar<sup>™</sup> is a wall mounted version of the award winning floor-to-ceiling SuperPole<sup>™</sup> with SuperBar<sup>™</sup>. It can be easily mounted to a wall much like a grab bar and combines all of the proven benefits of the SuperPole<sup>™</sup> with SuperBar's pivoting and locking dual rail technology. Designed for bathing, toileting, and other environments where falls can occur.

## Features and Benefits

## 5 Locking Positions.

Lift to unlock dual rail, keep rail up to move freely, let rail down to securely lock in place. Stand to shower, lean over to adjust shower controls or use when lowering into the tub. Easy to use ladder-like support for getting on / off the toilet.

### Safe Bathtub Entry.

Dual rail extends 20" / 51 cm from wall to provide a weight bearing surface for safe entry and exit from the bathtub.

### Easily Installs Like A Grab Bar.

Mounts securely to structural studs, with 3 contact points for added stength. Bottom grab bar is reversible for left or right side installation.

#### Stainless Steel Construction.

Unlike zinc or chrome plated grab bars that can chip, rust and stain tile over time, the stainless steel Dependa-Bar™ will never rust.





# **Technical Specifications:**

Dependa-Bar™	20" / 51 cm dual rail length x 1.25" / 3.175 cm diameter (ADA Compliant). Locks every $45^{\circ}$ (5 positions).
Rail Spacing	5" / 13 cm spacing between rails.
Rail Height	33" / 84 cm.
Lower Grab Bar	18" / 46 cm or 32" / 81 cm.
Wall Clearance	1.5" / 3.81 cm from wall (ADA Compliant).
Material	Stainless steel.
Finishes	White powder coat or knurled stainless steel.
Weight Capacity	325 lbs / 147kg.

