

**REPLACEMENT ACCESSORIES**

**SEAT AND BACKREST COVERS**

- 460900252 Cover-Set Comfort, White, 1/bx
- 460900312 Cover-Set Comfort, Blue, 1/bx
- 460900432 Cover-Set Comfort, Grey, 1/bx

**SPARE PARTS**


- 460900403 Charger for Bellavita, 1/bx
- 415252 Plug + Washer for Cover, White, 1/ea
- 415312 Plug + Washer for Cover, Blue, 1/ea
- 415432 Plug + Washer for Cover, Grey, 1/ea
- 460900213 Suction Cup with clip, 2.95" (75 mm) Grey, 1/ea
- 460900220 Set of Suction Cups with Clips, 4 x 2.95" (75 mm) Grey, 4/pk
- 460900600 Hand Control Including Storage Battery, 1/ea

**ACCESSORIES**

- 101385 Rotating & Transfer Aid, White, 1/cs
- 102385 Rotating Aid, White, 1/cs
- 410200252 Headrest, Padded, White, 1/ea
- 410200312 Headrest, Padded, Blue, 1/ea
- 410200432 Headrest, Padded, Grey, 1/ea
- 460900205 Side Flap, 1/ea



**bellavita**



**HIGHEST SEAT POSITION**



**LIGHTEST BATHLIFT AVAILABLE**



**RECLINING BACKREST LOWEST BATHLIFT AVAILABLE**

**SPECIFICATIONS:**

Overall Length (with backrest in a reclined position)	32.6"
Overall Width (seat incl. side flaps)	27.5"
Width of Seat	13.7"
Depth of seat	16.5"
Overall Length of Seat	22.4"
Sitting Height Min.	2.3"
Sitting Height Max.	18.8"
Height of Backrest	25.1"
Width of Backrest	12.5" - 14.1"
Backrest Reclines	50°
Width of Baseplate	12.2"
Length of Baseplate	22"
Load Capacity Max.	300 lbs.
Total Weight (excl. hand control)	20.5 lbs.
Shipping Weight	25.3 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

**bath safety**

To order call toll free: 877.224.0946

**drive**  
WWW.DRIVEMEDICAL.COM