



**ENVIRONMENT FRIENDLY PRODUCT**

Recycled materials were used to make this product.



**THIS ITEM IS AVAILABLE IN A 1-PACK!**

**K.D. Combination Plastic Transfer Bench/Commode**

- 12011KDC-2 Knock Down, Tool Free Removable Back and Arm, (Has replaced 12011C-2), 2/cs
- 12011KDC-1 Knock Down, Tool Free Removable Back and Arm, 1/cs
- 12011BRKT Replacement Commode Bracket for 12011C-2, 1/ea
- 12011SS Replacement Suction Cup for Transfer Bench, Small, 1pr/bx
- 12011SL Replacement Suction Cup for Transfer Bench, Large, 1pr/bx

**FEATURES:**

- Combines a transfer bench and commode into one product.
- New "A" frame construction provides additional stability.
- Durable blow-molded plastic bench and backrest.
- No exposed hardware to injure patient.
- Extra large, locking suction cups provide added safety.
- Comes complete with a 7.5 qt. commode bucket with metal handle, cover and splash shield.
- Limited Lifetime Warranty.



**K.D. Combination Padded Transfer Bench/Commode**



**THIS ITEM IS AVAILABLE IN A 1-PACK!**

- 12005KDC-2\* Combination Padded Transfer Bench/Commode, (Has replaced 12005C-2), 2/cs
- 12005KDC-1 Combination Padded Transfer Bench/Commode, (Has replaced 12005C-2), 1/cs
- 12011BRKT Replacement Commode Bracket for 12005C-2, 1/ea
- 12011SS Replacement Suction Cup for Transfer Bench Small, 1pr/bx
- 12011SL Replacement Suction Cup with Screws for Transfer Bench, Large, 1pr/bx

**FEATURES:**

- Combines a transfer bench and commode into one product.
- Comfortable cushioned seat and backrest.
- 1" Aluminum frame is lightweight, sturdy and corrosion resistant.
- Comes with 7.5 qt. commode bucket.
- Height adjusts in 1" increments.
- Back reverses without tools.
- Extra large suction cups provide added safety.
- Limited Lifetime Warranty.



**SUPPORT COLLAR REDUCES RATTLE!**  
**SPECIFICATIONS**

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12011KDC-1/2	28"	26.25"	17"	16.25"	21.75"	17.75"	10.3 lbs.	12.8/23.8 lbs.	400 lbs.
12005KDC-1/2	27.5"	24"	16.5"	16"	21.75"	17.75"	11.2 lbs.	13.7/25.4 lbs.	400 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.