



**SUPPORT COLLAR PREVENTS LEG MOVEMENT!**  
12203KD-4

**NEW PINCH FREE COVER ALLOWS FOR PUSH PINS TO BE DEPRESSED WITHOUT PINCHING FINGER. THE AUDIBLE "CLICK" WILL BE HEARD WHEN SECURE!**



**ENVIRONMENT FRIENDLY PRODUCT**  
Recycled materials were used to make this product.

**Deluxe, K.D. Aluminum, Bath Bench without Back** **Available 1 Pack**

- 12203KD-4 Knocked Down, Tool Free Assembly, (Has replaced 12106-4 and 12203-4), 4/cs
- RTL12203KDR Knocked Down, Tool Free Assembly, Retail, (Has replaced 12203KDR-1), 1/cs
- 12203KDRB-1 Knocked Down, Tool Free Assembly, Retail, Blue, 1/cs
- 12203KDRT-1 Knocked Down, Tool Free Assembly, Retail, Teal, 1/cs
- 12203KDRY-1 Knocked Down, Tool Free Assembly, Retail, Yellow, 1/cs
- 12002T Replacement Tips for 12203, 4/bx

**FEATURES:**

- Easy snap together assembly, Seat and Legs.
- Blow-molded bench provides comfort and strength.
- Drainage holes in bench reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- Height adjustable, angled legs with suction style tips provide additional stability.
- Support collar reduces rattle.
- Limited Lifetime Warranty.

**SPECIFICATIONS**

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12203KD-4	17"	20"	17"	12"	19.75"	15.75"	4.5 lbs.	21 lbs.	400 lbs.
RTL12203KDR	17"	20"	17"	12"	19.75"	15.75"	5 lbs.	7 lbs.	400 lbs.
12203KDRB, T, Y-1	17"	20"	17"	12"	19.75"	15.75"	5 lbs.	7 lbs.	400 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.