



**TOOL FREE  
REMOVABLE  
BACK, LEGS  
& ARM!**

**BARIATRIC**  
**500 LB**  
WEIGHT CAPACITY



**ENVIRONMENT FRIENDLY PRODUCT**  
Recycled materials were used  
to make this product.

Available  
**1 Pack**

12021C



**EURO-STYLE  
CLIP ENSURES  
A SECURE FIT!**



12025KD-2



**MUST PULL "TAB" ON  
SUCTION CUP TO RELEASE  
FROM THE TUB.**

### K. D. Bariatric Transfer Bench

- 12025KD-2 Tool Free Assembly, (Has replaced 12025-2), 2/cs
- 12025KD-1 (Has replaced 12025-1 and 12009-1), 1/cs
- 12011SS Replacement Suction Cup for Transfer Bench, Small, 1pr/bx
- 12011SL Replacement Suction Cup for Transfer Bench, Large, 1pr/bx

#### FEATURES:

- 500 lb. Weight Capacity.
- "A" frame construction provides additional stability.
- Durable blow-molded plastic bench and backrest.
- Aluminum frame is lightweight, sturdy and corrosion resistant.
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs.
- All transfer benches are reversible to accommodate any bathroom.
- No exposed hardware to injure patient.
- Extra large, locking suction cups provide added safety.
- Tool-Free Assembly requires connecting back, legs and arm.
- Limited Lifetime Warranty.



12009-2

**EXTRA MIDDLE LEGS  
FOR ADDED STRENGTH.**

Available  
**1 Pack**

**BARIATRIC**  
**600 LB**  
WEIGHT CAPACITY

### Deluxe Bariatric Transfer Bench with Cross Frame Brace

- 12009-2 Fully Assembled, 2/cs
- 12011SS Replacement Suction Cups for 12009, 6/bx

#### FEATURES:

- 600 lb. Weight Capacity.
- Durable, blow-molded plastic bench and backrest.
- Aluminum frame is lightweight, sturdy and corrosion resistant.
- Extra middle legs for added strength.
- Limited Lifetime Warranty.

#### SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12025KD-1/2	33"	26.5"	25"	16.5"	23.25"	18.25"	12 lbs.	14/27 lbs.	500 lbs.
12009-2	32"	27"	22.25"	16"	21"	16"	14 lbs.	32 lbs.	600 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.