



SIDE POUCHES MAKE IT EASY TO CARRY ITEMS WITH YOU!

COMES WITH GLIDE SKIS!

300 LB
WEIGHT CAPACITY

EASY FOLDING DESIGN USES ELASTIC CORDS!



RTL10263KDR COMES WITH 5" WHEELS!

Deluxe Folding Walker

RTL10263KDR Deluxe, Folding Travel Walker with 5" Wheels and Fold Up Legs, (Has replaced 10263KDR), 1/cs


Available **1 Pack**

FEATURES:

- Easily folds tool free to 1/2 the size of a traditional walker and can be stored in a "Ready Set Go" pouch provided with each walker.
- Trigger release feature allows the user to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity. Trigger release can be reached with hands of all sizes.
- Comes with 5" front wheels and rear glide skis that are ideal for most surfaces.
- Two side storage pouches with large pockets for personal items are held in place by Velcro tabs to both side frames.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight.
- Deluxe soft ribbed contoured hand grips.
- Limited Lifetime Warranty.



LEGS CONVENIENTLY FOLD UP FOR EASY STORAGE!




WALKER CAN BE FOLDED FLAT FOR STORAGE & TRANSPORT!



LEGS ARE ATTACHED WITH STRONG ELASTIC CORDS SO LEGS CAN BE FOLDED!



COMES WITH READY SET GO BAG!



PIN MUST GO INTO LOCK COLLAR BEFORE USE.

SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
RTL10263KDR	17"	22"	26.5"	20.5"	37"	29"	8 lbs.	10 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.