





12487





RTL12350KDR Triangular Seat, Knock Down,

Tool Free Assembly, Retail, (Has replaced 12350KDR-1), 1/cs

RTL12351KDR Rectangular Seat, Knock Down,

Tool Free Assembly, Retail, (Has replaced 12351KDR-1), 1/cs

10389G-1 Replacement Tips for 12350, 12351, 4/bx

FEATURES:

- Attractive high-quality teak is renowned for its durability and weathering properties.
- Ideal for use in shower stall.
- Teak will naturally take on a handsome silver grey patina after several months of exposure to water.
- Ideal for a woman's grooming stool.
- Lightweight, durable, adjustable legs.
- Large rubber tips reduce slipping.
- Limited Lifetime Warranty.





12486 Without Back, Retail, 1/cs 12487 With Back, Retail, 1/cs

12486T Replacement Tips for 12486, 12487, 4/bx

Pack

FEATURES:

- Conveniently folds flat.
- Blow-molded bench provides comfort and strength.
- Drainage holes in seat reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- Angled legs provide additional stability.
- Limited Lifetime Warranty.



CONVENIENTLY FOLDS FLAT!

SPECIFICATIONS									
Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
RTL12350KDR	22"	22.5"	22"	16"	19"	13"	9 lbs.	11 lbs.	300 lbs.
RTL12351KDR	17"	20"	17"	10"	18"	13"	7.5 lbs.	9 lbs.	300 lbs.
12486	17.5"	20"	20"	12"	18"	16"	4.5 lbs.	6.6 lbs.	300 lbs.
12487	17.5"	20"	20"	12"	18"	16"	5 lbs.	7 lbs.	300 lbs.

^{*}Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.

12486