



Economy Knee Walker

780 Adult, Knee Walker,

300 lb. Weight Capacity, 1/cs

780 BASKET Basket for 780 Knee Walkers, 1/ea 9503H78003 Seat for 780, Knee Walkers, 1/ea 9503H78007 Brake for 780, Knee Walkers, 1/ea

FEATURES:

- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot.
- Deluxe lever brake ensures safety.
- Tool free height adjustment.
- 8" casters are ideal for indoor/outdoor use.
- 300 lb. Weight Capacity.
- Limited Lifetime Warranty.
 Lifetime Warranty on brake cable.





780 BASKET (SOLD SEPARATELY FOR USE WITH 780 DRIVE KNEE WALKERS)



CONVENIENTLY FOLDS FOR STORAGE AND TRANSPORT



LEVER LOCK CAN BE ADJUSTED FOR EITHER RIGHT OR LEFT HANDED USE



LEG HEIGHT ADJUSTS TO FIT THE HEIGHT OF THE USER



TOOL FREE, AND EASY ADJUSTMENTS FOR RIGHT OR LEFT LEG SUPPORT



HANDLE ADJUSTS TO FIT THE HEIGHT OF THE USER

SPECIFICATIONS

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Depth	Pad Height (highest)	Pad Height (lowest)	Pad Width	Pad Depth	Product Weight	Carton Shipping Weight	Weight Limit
780	41 "	33"	14.5"	29"	22"	18"	7.5"	14.5"	16 lbs.	19 lbs.	300 lbs.

^{*}Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.