

# LIFESTYLE ESSENTIALS



## MOLDED STOCKING AID

RTL2010

3/cs

### Features:

- Makes it possible to put on stockings without bending.
- Curved cut outs on the sides hold the sock while it is being pulled up the leg.
- Stocking aid holds the sock open and allows the user to easily pull the sock or stocking up with minimal effort.
- Limited 6 Month Warranty.



## POLYETHYLENE STOCKING AID

RTL2012

3/cs

### Features:

- Soft plastic stocking aid that can be easily folded or rolled to fit inside a sock.
- Curved cut outs on the sides hold the sock while it is being pulled up the leg.
- A long webbing strap allows the user to pull the strap without bending over.
- Latex free.
- Limited 6 Month Warranty.



contact your Sales Representative, or visit [www.drivemedical.com](http://www.drivemedical.com)