Βάτη Σάγετν

oline



NEW PINCH FREE COVER ALLOWS FOR PUSH PINS TO BE DEPRESSED WITHOUT PINCHING FINGER. THE AUDIBLE "CLICK" WILL BE HEARD WHEN SECURE!


Available K.D. Padded Transfer Bench

12005KD-2	Knock Down Padded, Tool Free Assembly Back,
	Legs and Arm, (Has replaced 12005-2), 2/cs
12005KD-1	Knock Down Padded, Tool Free Assembly Back,
	Legs and Arm, (Has replaced 12005-1), 1/cs
12005KDR-1	Knock Down Padded, Tool Free Assembly
	Back, Legs and Arm, Retail, 1/cs
12011SS	Replacement Suction Cup for Transfer Bench,
	Small, 1pr/bx
12011SL	Replacement Suction Cup with Screws for
	Transfer Bench, Large, 1pr/bx
12005TP	Replacement Transfer Pad, 1/ea

FEATURES:

- Comfortable cushioned seat and backrest.
- "A" frame construction provides additional stability.
- Height adjusts in ¹/2" increments with unique "Dual Column" extension legs.
- All transfer benches are reversible to accommodate any bathroom.
- Aluminum frame is lightweight, sturdy and corrosion resistant.
- No exposed hardware to injure patient.
- Extra large suction cups provide added safety.
- Tool-Free Assembly requires connecting back, legs and arm.
- Limited Lifetime Warranty.







SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12005KD-2	27.5 "	24″	16.5 "	16″	21.75"	17.75″	11.75 lbs.	26.5 lbs.	400 lbs.
12005KD-1`	27.5″	24″	16.5 "	16 "	21.75″	17.75″	11.75 lbs.	14 lbs.	400 lbs.
12005KDR-1	27.5 "	24 "	16.5 "	16″	21.75"	17.75″	11.75 lbs.	16 lbs.	400 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.