



Available K.D. Plastic Transfer Bench

| 12011KD-2 | Knock Down Tool Free Back, Legs and Arm, 2/cs (Has replaced 12011-2 and 12020-2) |
|-------------|---|
| 12011KD-1* | Knock Down Tool Free Back, Legs and Arm, 1/cs (Has replaced 12011-1) |
| RTL12011KDR | Knock Down, Retail, Tool Free Back, Legs and Arm, (Has replaced 12011KDR-1), 1/cs |
| 12011SS | Replacement Suction Cup for Transfer Bench, Small, 1pr/bx |
| 12011SL | Replacement Suction Cup with Screws for Transfer Bench, Large, 1pr/bx |
| 12011-A | Replacement Arm for 12011KD and 12011 1/ea |

- "A" frame construction provides additional stability.
- Durable blow-molded plastic bench
- Height adjusts in 1/2'' increments with unique "Dual Column" extension legs.
- Aluminum frame is lightweight, sturdy and corrosion resistant.
- All transfer benches are reversible to accommodate any bathroom.
- No exposed hardware to injure patient.
- Extra large, locking suction cups provide
- Tool-Free Assembly requires connecting back,
- Limited Lifetime Warranty.

$(\mathbf{1})$





SPECIFICATIONS

| Item # | Width (outside legs) | Width Seat | Depth (outside legs) | Depth Seat | Seat Height (highest) | Seat Height (lowest) | Product Weight | Carton Shipping Weight | Weight Limit |
|-------------|----------------------------|---------------|----------------------------|-----------------|-----------------------------|----------------------------|-------------------|------------------------------|-----------------|
| 12011KD-2 | 33 " | 26 " | 24 " | 18.5 " - 19.5 " | 21.75 " | 17.75″ | 10.5 lbs. | 24 lbs. | 400 lbs. |
| 12011KD-1 | 33 " | 26 " | 24″ | 18.5 " - 19.5 " | 21.75″ | 17.75″ | 10.5 lbs. | 13 lbs. | 400 lbs. |
| RTL12011KDR | 33 " | 26 " | 24 " | 18.5" - 19.5" | 21.75" | 17.75″ | 10.5 lbs. | 14 lbs. | 400 lbs. |

16

BATH SAFETY



Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.