



Suction Cup Grab Bars




THIS ITEM IS AVAILABLE IN A 1-PACK!

- 13082-1 12" Grab Bar, 1/cs
- 13082-3 12" Grab Bar, Retail, 3/cs
- 13084-1 19 3/4" Dual Grab Bar, Rotating, 1/cs
- 13084-3 19 3/4" Dual Grab Bar, Rotating, Retail, 3/cs

FEATURES:

- Installs and is removable without tools or professional installation.
- Will not damage property.
- Large suction cups provide an extremely strong hold.
- Release levers make installing and removing the suction cup grab bars quick and easy.
- Color indicator shows "safe secure" suction or "re-affix" suction.
- Limited Lifetime Warranty.



THE SUCTION CUP GRAB BAR ROTATES TO ANY ANGLE!



RED/GREEN COLOR INDICATORS SHOW "SAFE SECURE" SUCTION OR "RE-AFFIX" SUCTION.

CAUTION: THE SUCTION CUP GRAB BARS ARE TO BE USED AS BALANCE ASSIST ONLY. THEY ARE NOT FOR BODY WEIGHT LEVERAGE.

SPECIFICATIONS

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

| Item # | Length | Product Weight | Shipping Weight |
|---------|--------|----------------|-----------------|
| 13082-1 | 12" | .4 lbs. | .9 lbs |
| 13082-3 | 12" | .4 lbs | 1.7 lbs |
| 13084-1 | 19.75" | 1.0 lbs | 1.5 lbs |
| 13084-3 | 19.75" | 1.0 lbs | 3.5 lbs |

bath safety

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM