

AQUATEC®

INVACARE® CANADA
Aquatec® Ocean E-VIP



Yes, you can.®

The Aquatec Ocean E-VIP electric tilting commode is the newest in a long line of innovative Ocean Commode products. Designed with the healthcare provider in mind, the E-VIP reduces back strain, protects the personal hygiene of the attendant, and limits the need for physically lifting the patient, all while ensuring that the highest level of safety and comfort are provided for the patient at all times.

Features

- Electrically adjustable tilting-angle for increased safety
- High weight capacity (330 lbs.), stainless steel frame
- Tilting-angle can be adjusted between -5° and 35°
- Ergonomic seating is accompanied by a wide range of accessories
- Seating height can be electrically adjusted between 20" - 41"
- Tremendous stability even in most elevated positions
- Large space between side tubes for use as toilet commode
- Complete chair can be machine cleaned at up to 85°C
- Battery capacity for 40 complete lifting cycles with full load
- Battery can be charged separately
- Compatible with most Ocean/ Ocean VIP accessories*

*Except reclining backrest and anti-tippers

Model no. OCEANEVIP, Aquatec Ocean E- VIP Commode

Specifications

Seat Depth	18"
Seat Width	19"
Seat Height	20 - 41"
Overall Width	29"
Overall Length	40"
Overall Height	42 - 64"
Distance from Seat to Footrest	13 1/2 - 18"
Front Seat to Floor	16 1/2 - 38 1/2"
Back Seat to Floor	18 - 40"
Width between armrests	18"
Distance between rear wheels	25"
Tilt angle forward	-5°
Tilt angle backward	35°
Weight capacity (lbs.)	330
Weight (lbs.)	86

Aquatec Ocean shower commodes are covered under a two year limited warranty against manufacturing defects. Manufacturer's warranty does not cover normal wear.



Invacare Canada LP
www.invacare.ca

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2012
Invacare Canada

All rights reserved.
Trademarks are identified by the symbols TM and ®. All trademarks are owned or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.

Form No. 08-721C
Rev. 03/12

