









A quatec Shower Commodes offer superior functionality, maximum adjustability and an enhanced quality of life. Aquatec commodes feature an innovative, practical design that is ideal for self-propelling clients as well as homecare and nursing facility applications.

For safe and proper use:

- Ensure that the brakes are locked in parked position when in use
- Follow Aquatec's care instructions to keep your commode clean
- Can be used as aids for all bathroom needs
- For indoor use only



Aquatec, your first choice for safety and reliability!

Aquatec® Shower Commodes



Ocean Standard/Ocean SP

Model no. 7.10.900K Standard
Model no. 7.12.900K SP (Self Propelled)

The multi-functional Ocean Standard and Ocean SP commodes are perfect for institutional or home use. With maximum adjustability, both can be used as a commode chair, shower chair or as a bedside commode. The Ocean SP also features 24" rear wheels making it ideal for self-propelling.

Features

- Tool-less seat height adjustment range of 5" to fit over most toilets
- Removable height-adjustable swing-away footrests with adjustable heel loops
- Easy to clean smooth surfaces, detachable seat and tensionadjustable, machine-washable backrest
- Four, 5" casters with wheel & swivel locking mechanism
- Rust-proof frame ideal for pool use
- Commode pan with lid included
- Standard with 24" rear wheels with SP model

Specifications

	Standard	SP
Model	7.10.900K	7.12.900K
Seat Depth	18"	18"
Seat Width	19"	19"
Seat Height	19" - 24"	19" - 24"
Overall Width	22"	26 ½"
Overall Length	36 ½ "	41 ½"
Seat to Floor	19" - 24"	19" - 24"
Height Under Seat with Pan	13 ½" - 18 ½"	13 ½" - 18 ½"
Overall Height	38" - 43"	38" - 43"
Width Between Armrests	18"	18"
Distance Between Rear Wheels	N/A	21"
Stationary Posterior Tilt	5°	5°
Anterior Tilt	0°	0°
Weight Capacity (lbs.)	285	285
Product Weight (lbs.)	32	43

Aquatec® Shower Commodes





The XL offers 3" of additional space between the armrests

Ocean XL

Model no: 7.13.900K

The multi-functional Ocean XL has all of the great features of the Standard Ocean model with an increased weight capacity of 330 lbs. and 3" of additional space between the armrests for added comfort.

Features

- Sturdy stainless steel frame with increased weight capacity
- Tool-less seat height adjustment range of 5" to fit over most toilets
- Removable height-adjustable swing-away footrests with adjustable heel loops
- Easy to clean smooth surfaces, detachable seat and tensionadjustable, machine-washable backrest
- Four, 5" casters with wheel & swivel locking mechanism
- Rust-proof frame ideal for pool use
- Commode pan with lid included
- Available with 24" self-propelling wheels (part no. 16379)

Specifications

	Ocean XL
Model	7.13.900K
Seat Depth	18"
Seat Width	19"
Seat Height	19" - 24"
Overall Width	25, 26 ½"
Overall Length	37, 42"
Seat to Floor	19" - 24"
Height Under Seat with Pan	13 ½" - 18 ½"
Overall Height	38" - 43"
Width Between Armrests	21"
Distance Between Rear Wheels	21"
Stationary Posterior Tilt	5°
Anterior Tilt	O°
Weight Capacity (lbs.)	330
Product Weight (lbs.)	32

Aquatec® Shower Commodes



Ocean VIP Tilt and Soft Seat

Model no: 1470793KIT

The Ocean VIP Tilt shower commode offers maximum versatility with a stainless steel frame, various height adjustments and an advanced tilting function.

Features

- 35° posterior tilt
- Capability for multiple adjustments to accommodate the growth of a pediatric client and the changing needs of an adult user
- Easy to clean smooth surfaces, detachable seat and tensionadjustable, machine-washable backrest
- Structured fold-away, angle adjustable footplates with heel loops
- Height and depth adjustable headrest
- Armrest can be mounted in 2 different heights: 8" and 9 ½" from the seat
- Soft seat (with hygiene recess), commode pan and pan holder included

Specifications

	Ocean VIP
Model	1470793KIT
Seat Depth	18"
Seat Width	19 ½"
Seat Height	19 ½" - 23 ½"
Overall Width	22 ½"
Overall Length	41 ½"
Seat to Floor	19 ½" - 23 ½"
Height Under Seat with Pan	13 ½" - 18 ½"
Overall Height	43" - 47"
Width Between Armrests	18"
Distance Between Rear Wheels	28" - 34"
Stationary Posterior Tilt	5°
Anterior Tilt	0°
Weight Capacity (lbs.)	330
Product Weight (lbs.)	47

Aquatec® Shower Commode Accessories

Accessories



Sanitary Pan
With grip & lid.
Model no. 10230



Safety Bar
For additional security. Easy assembly.
Padded grip.
Model No. 10167-10
Available for Ocean Standard, SP & VIP only



Calf Support
Left or right side. Adjustable in height
and angle. Upholstered pad.
Model no. 16615
Available for Ocean Standard, SP & XL only



Frontal Splash Guard
Easily fitted. For soft seat with
hygiene recess.
Model no. 16511



Armrest Lock
Locks the armrest in the lower position.
Model no. 16597
Available for Ocean Standard, SP, XL only
Model no. 1470077
Available for Ocean VIP only



Soft Seat
With hygiene recess
Model no. 16341
Available with Ocean Standard, SP, XL only (shown)
Model no.1470075
Available for Ocean VIP only
With toileting recess
Model no. 16342
Available for Ocean Standard, SP, XL only

Aquatec[®] Commode Accessories & Warranty

Accessories



Amputation Support

Left or right side. Adjustable in height and angle. Upholstered pad. Model no. 16616 Available for Ocean Standard, SP, XL only



Insert for Soft Seat

For soft seat with hygiene recess. Model no. 16343



Anti-Tippers

For additional safety.
Model no. 1470833

Warranty

Aquatec Ocean shower commodes are covered under warranty for a period of two years.





Invacare Canada LP

www.invacare.ca

Canada

570 Matheson Blvd. E. Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

All rights reserved.
Trademarks are identified by the symbols ™, ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Specifications are subject to change without notification.

© 2009

