

MAX TRUEFITT BACK



Max Series Backs **STYLE: Max TrueFitt Back**
Stock Sizes – Same Day Shipment
ADP CODE: SEBFM0050

Part Number	Size W x H	Part Number	Size W x H
MTFB1616	16" x 16"	MTFB1820	18" x 20"
MTFB1620	16" x 20"	MTFB2016	20" x 16"
MTFB1816	18" x 16"	MTFB2020	20" x 20"


Non-Stock Sizes – Shipment within 3 business days
ADP CODE: SEBFM0055

Part Number	Size W x H	Part Number	Size W x H
MTFB1414	14" x 14"	MTFB1818	18" x 18"
MTFB1618	16" x 18"	MTFB1916	19" x 16"
MTFB1716	17" x 16"	MTFB1920	19" x 20"
MTFB1720	17" x 20"	MTFB2018	20" x 18"

Custom Sizes - Shipment within 7 business days
ADP CODE: SEBFM0065

Part Number	Size W x H
MTFBCUS	Custom

Options		
Part Number	Size W x H	Description
MTFBV	2" x 2"	Visco Foam
MTFB3	3" x 2"	Poly Foam




ADP CODE: SEBFM0150 - Max Back

(The upcharge of \$100 in addition to the ADP codes for the Ultra, Air, or TruFitt Backs)

Product Weight: 5.5 lbs (18"W x 16"H)



FUTURE
MOBILITY HEALTHCARE



Lateral Dimensions		Max Truefitt
MAX TRUEFITT, Heights 14, 16, 18, 20		
Chair Width	DIM B (")	DIM A (")
14"	10.5	7.6
16"	12.5	7.6
18"	14.5	7.6
20"	16.5	7.6
22"	18.5	7.6



- ✓ Fully adjustable for height, depth and angle of recline
- ✓ Hardware fits 3/4", 7/8" and 1" tubing and can be installed within a few minutes with the right tools
- ✓ Quick release lower latch levers allows the backrest to be easily removed.