



# ECLIPSE

by **pdg**  
MOBILITY TECHNOLOGIES

# bariatric

manual wheelchair



- For people of significant size
- Delivers style and performance
- Forward wheelbase for excellent maneuverability
- Transportable
- Narrow overall width



WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

[www.pdgmobility.com](http://www.pdgmobility.com)

# ECLIPSE

THE RIGHT MANUAL BARIATRIC WHEELCHAIR FOR YOUR NEXT CLINICAL APPLICATION.



600 lb capacity;  
1000 lb upgrade option

Unique wheelbase design ensures most of the weight is carried on the rear wheels, allowing front casters to pivot and roll with ease

Easily transportable with removeable back, rear wheels, armrests and front rigging

Rigid frame provides high strength and stiffness

Forward wheelbase

Blue Shimmer



Tracer Red



Black Shimmer



Burnt Orange



Forest Green



Misty Purple



TREND-SETTING FRAME COLORS. YOU CAN ORDER CUSTOM COLORS TOO. CONTACT US FOR DETAILS.



Transportable.



Forward wheelbase.



Adjustable seat depth and back angle.



Narrow overall width.



Solar cushions specially designed for bariatric applications.

Seat width (adjustable) 20" - 48"  
Seat depth (adjustable) 16" - 24"  
Seat height (adjustable) 14" - 20"  
Back angle (adjustable) 85° - 100°  
Recliner Reclines to 130°  
Back height 25" and 30"  
Overall width Seat width plus 7"

Overall length (no front rigging) 31.5" - 34.5"  
Weight capacity 600 lbs or 1000 lbs upgrade  
Front rigging 60° Swing-away, elevating, center mount  
Chair weight 35 lbs (without backposts, armrests, front rigging, rear wheels, anti-tippers)  
Warranty PDG's exclusive TGIF warranty (Thank Goodness It Fits)  
Lifetime warranty on the frame.

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product. Canadian Patent # 2,198,643 US Patent # 5,957,474

**USA**  
PDG Product Design Group  
PMB # 173 - 1124 Fir Avenue  
Blaine, WA  
98230, USA

**CANADA**  
PDG Product Design Group  
Unit 102 - 366 East Kent Avenue South  
Vancouver BC  
V5X 4N6

T: +1 604 323 9220  
F: +1 604 323 9097  
Toll Free (in N. America): 888-858-4422  
Toll Free Fax (in N. America): 877-323-9093  
www.pdgmobility.com