



# FUZET20

by **pdg**  
MOBILITY TECHNOLOGIES

20° of

*tilt*



Forward pivot point provides 20° of tilt with no knee rise

- Tilt for positioning and improved posture
- Foot and hand propel in tilted position
- Maneuverable short wheelbase
- Width adjustable
- Depth adjustable
- Transit ready option


**pdg**  
MOBILITY TECHNOLOGIES

WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

[www.pdgmobility.com](http://www.pdgmobility.com)

# FUZET20

THE RIGHT MANUAL TILT-IN-SPACE WHEELCHAIR FOR YOUR NEXT CLINICAL APPLICATION.



TREND-SETTING FRAME COLORS. YOU CAN ORDER CUSTOM COLORS TOO. CONTACT US FOR DETAILS.




Adjustable back angle



Available with self tilt option



Optional transit tie-downs



Adjustable axle and seat height



Adjustable chair width

Seat width (adjustable) 13" – 20"  
 Seat depth (adjustable) 15" – 20"  
 Seat height (adjustable) 13" – 20"  
 Back angle (adjustable) 90° – 120°  
 Back height 20", 25"  
 Dynamic tilt -5° to 20°

### Overall width - T20 Standard

- for 12" & 16" rear wheels, use seat width plus 9.5"  
 - for 20", 22" & 24" rear wheels, use seat width plus 10.75"

### Overall width - T20 Narrow

- for 12" & 16" rear wheels, use seat width plus 7.5"  
 - for 20", 22" & 24" rear wheels, use seat width plus 8.75"

Overall Length (no front rigging) 28"  
 Weight Capacity 250 lbs  
 Chair Weight 35 lbs (without front rigging, wheels, anti-tippers, armrests)  
 Transit Options Available ANSI/RESNA WC-19  
 Warranty PDG's exclusive TGIF warranty (Thank Goodness It Fits) Lifetime guarantee on the frame.

**USA**  
 PDG Product Design Group  
 PMB # 173 - 1124 Fir Avenue  
 Blaine, WA  
 98230, USA

**CANADA**  
 PDG Product Design Group  
 Unit 102 - 366 East Kent Avenue South  
 Vancouver BC  
 V5X 4N6

T: +1 604 323 9220  
 F: +1 604 323 9097  
 Toll Free (in N. America): 888-858-4422  
 Toll Free Fax (in N. America): 877-323-9093  
[www.pdgmobility.com](http://www.pdgmobility.com)