



FUZET50JR

by **pdg**
MOBILITY TECHNOLOGIES

effortless **tilt**

50° of effortless tilt,
zero center-of-gravity shift



- Small Footprint
- ISO tested
- WC-19 Crash Tested
- Transit Ready option
- Grows with ease


pdg
MOBILITY TECHNOLOGIES

WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

www.pdgmobility.com

FUZET 50JR

JUNIOR TILT-IN-SPACE AT ITS BEST!



TREND-SETTING FRAME COLORS. YOU CAN ORDER CUSTOM COLORS TOO. CONTACT US FOR DETAILS.



Adjustable seat height



Adjustable axle position



Adjustable seat depth and back rest



Transit tie downs



Adjustable seat width

Seat width (adjustable) 12"–18"
 Seat depth (adjustable) 13"–18"
 Seat height (adjustable) 13"–20"
 Back angle (adjustable) 90°–120°
 Back height 20" and 25"
 Dynamic tilt 50°

Overall width
 - for 12" & 16" rear wheels, use seat width plus 7.5"
 - for 20", 22" & 24" rear wheels, use seat width plus 9"

Overall Length (no front rigging) 24"
 Weight Capacity 180 lbs
 Chair Weight 31 lbs (without front rigging, wheels, anti-tippers, armrests)
 Transit Options Available ANSI/RESNA WC-19
 Warranty PDG's exclusive TGIF warranty (Thank Goodness It Fits) Lifetime guarantee on the frame.

USA
 PDG Product Design Group
 PMB # 173 - 1124 Fir Avenue
 Blaine, WA
 98230, USA

CANADA
 PDG Product Design Group
 Unit 102 - 366 East Kent Avenue South
 Vancouver BC
 V5X 4N6

T: +1 604 323 9220
 F: +1 604 323 9097
 Toll Free (in N. America): 888-858-4422
 Toll Free Fax (in N. America): 877-323-9093
www.pdgmobility.com