



STELLAR

by **pdg**
MOBILITY TECHNOLOGIES

manual *tilt*



45° of tilt
for pressure
relief and positioning

- Ideal for hand propulsion with accessible rear wheels
- Knees remain low in tilt - wheel under desks, into vans and access other environments with ease

- Short wheelbase provides excellent maneuverability
- Features adjustable seat height, seat depth and back angle
- Available in Heavy Duty (HD) extra-wide sizes



Stellar HD

pdg
MOBILITY TECHNOLOGIES

WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

www.pdgmobility.com

STELLAR

THE RIGHT MANUAL TILT-IN-SPACE WHEELCHAIR FOR YOUR NEXT CLINICAL APPLICATION.



Blue Shimmer



Tracer Red



Black Shimmer



Burnt Orange



Forest Green



Misty Purple



TREND-SETTING FRAME COLORS. YOU CAN ORDER CUSTOM COLORS TOO. CONTACT US FOR DETAILS.



Multi-adjustable



Adjustable push bar option



Easily adjustable axle position



Variety of caster and wheel sizes



HD option

Seat width, Standard 14" – 22"
Seat width, HD 23" – 32"
Seat depth (adjustable) 16" – 22"
Seat height (adjustable) 14" – 20"
Back angle (adjustable) 90° – 120°
Back height (adjustable) 18" – 26"
Dynamic tilt 45°

Overall width
 - for 12" & 16" rear wheels, use seat width plus 9.5"
 - for 20" & 24" rear wheels, use seat width plus 10.5"

Overall length (no front rigging) 30" – 37"
Weight capacity 250 lbs, HD to 450 lbs (Custom to 600 lbs)
Chair weight 40 lbs (HD 50lbs frame weight)
Transit option available
Warranty ANS/RESNA WC-19 (Standard frame only)
 PDG's exclusive TGIF warranty (Thank Goodness It Fits)
 Lifetime warranty on the frame.

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product. US Patent # 6 447064 B1

USA
 PDG Product Design Group
 PMB # 173 - 1124 Fir Avenue
 Blaine, WA
 98230, USA

CANADA
 PDG Product Design Group
 Unit 103 - 318 East Kent Avenue South
 Vancouver BC
 V5X 4N6

T: +1 604 323 9220
 F: +1 604 323 9097
 Toll Free (in N. America): 888-858-4422
 Toll Free Fax (in N. America): 877-323-9093
www.pdgmobility.com