

The logo for QUICKIE, featuring a stylized orange swoosh above the word "QUICKIE" in white, bold, sans-serif capital letters.

X2

MID-WHEEL DRIVE

The New Mid-Wheel Drive Power Wheelchair



MADE IN CANADA 



is the new addition to the Canadian-made power wheelchair family. This new mid-wheel drive wheelchair suggests a unique offering for a user sharing their time indoors and outdoors. The low seat to floor height and the short turning radius will facilitate day to day activities at home, at play and at work. The Quickie® X2™ will clear outdoor obstacles for a smooth and dependable ride.



Mid-Wheel Drive Featuring:

- GC3 Techonology
- Short Turning Radius
- PG Drives Technology
- Fully Adjustable Rehab Seat
- Low Seat to Floor Height





Mid-Wheel Drive 2 Pole Power Wheelchair

The Quickie® X2™ has a unsurpassed maneuverability which revolutionizes the standards in power wheelchairs.

GC³ TECHNOLOGY™



1. Gravity Control

Stability and traction without necessary adjustments to compensate for various body types relative to complex or dynamic seating.



2. Ground Compensation

While riding over uneven surfaces, the hard movements and vibrations are controlled and not absorbed by the seat or user.



3. Gravity compensation

The unique suspension will distribute the weight over the entire footprint and assure an anterior and posterior stability.

SEATING

Fully Adjustable Rehab Seat

- Adjustable Seat width 15" to 20"
- Adjustable Seat Depth 15" to 19"
- Pre-Set back angle: 86° to 122°- increments of 4°
- Eclipse seat optional



Complement your Quickie® X2™ with JAY J3 seating

FEATURES

Weight Capacity	300 lbs Standard
Front Caster	6" Caster
Drive Wheel	12.5" Cast wheel pneumatic with or without insert
Transit	4-Point wheelchair transport included
Electronics	Standard: VR2 70 Amp Options: VR2-90 Amp, R-Net EL 90 Amp expandable, LED joystick
Motor	2 Pole
Seat to floor height	15 3/4"
Max Speed	5.6 MPH
Overall Base Width	24"
Overall Base Length	43.5"
Seating	Fully adjustable seat: Eclipse
Warranty	5 years for frame, please see owner's manual for all warranty information
Colours available	Standard: Black Options: Pearl Red and Pearl Blue



www.QuickieXPower.ca

Quickie® X2™ REV.C



Customer Service: 1.905.660.2459 / 1.800.263.3390
or visit www.SunriseMedical.ca

improving people's lives