







EASILY GET OUT OF BED

EASILY PIVOIS DOWN

Features & Benefits

- Bed Rail/Hand Rail Adjustable length allows for use as both a bed rail to prevent falling out of bed or as a hand rail to asist in getting out of bed.
- Adjustable Bed rail easily adjusts from 26" to 34" to 42" after it has been placed on the bed.
- Collapsible Rail folds down to side of bed to allow user to get out of bed or provide space for making bed.
- Reversible Can be placed on either side of the bed
- Dual Safety Strap Secures bed rail to any home or hospital bed. Accommodates up to King size beds.

Specifications

Width of Rail: Locks in place at 26", 34" & 42" after installed on bed Height of Rail: 23" high from the base of mattress

Gap Between Rails: 3"

Bed Type: Attaches to any home or hospital bed with included Safety Strap

Weight of Product: 14 lbs Weight Capacity: 400 lbs

Assembly: Installs in minutes with 4 bolts and Allen Wrench

(Included)

Package Dimensions: 28" x 21" x 1.5"

Case of 4 Dimensions: 28" x 21" x 6"; 60 lbs

How Do We Compare	EZ Adjust Bed Rail	Other Bed Rails
Adjustable Length Before and After Installation		&
Pivots Out of the Way	Ø	×
Versatility (Bed Rail and Bed Handle)	Ø	×
Accommodates All Bed Types	©	×
Durable Powder Coat Finish	Ø	⊘
Prevents Falling from the Bed	Ø	⊘

Easy Assembly (No Additional Tools Required)



Attach L Shaped Tubes



Slide On Strap



Attach Rear Crossbar



Place On Bed