

KID KART

family album



Kid Kart cares
about Kids

KIDKART Xpress

Stage One



Height-adjustable push handle

Standard Sunshade

Ability to position shell forward or rearward facing

Light weight, high strength frame is easy to fold and lift for transport

Multiple quick release rear wheel and front caster options make travel easy over any terrain

Convenient storage compartment carries up to 15 lbs.




Portable

Folds with or without the seating system.


Kid Kart Xpress

Configuration shown with shell facing caregiver, Purple Frame, Jazz Fabric Package and Navy Seating.



Transit Option

Special Features



Height-Adjustable Booster Base

*Height-adjustable floor
sitter with locking casters
collapses flat for storage.*

Raincape Option

*Water resistant fabric-
boot and hood can be
used separately.*



Frame Color Options



PEARL PINK



PURPLE




TEAL

Fabric Design Options



CLASSIC



JAZZ



ORIGINAL





KIDKART TLC

Cool Features & Benefits



Comfortable, foam-padded, height-adjustable handlebar for varying caregiver needs

Sunshade

Ability to position shell forward or rearward facing

Removable storage pack

Quick-release rear wheels



Fixed height feeder base with locking casters.

Frame Color Options

Fabric Design Options



DUSTY PLUM



BLACK



ROYAL BLUE and FUSCIA with SILVER TRIM



FOREST GREEN and PURPLE with SILVER TRIM



PURPLE and TEAL with FUSCIA TRIM



KID KART Custom Seating System

Tri-wing head support (shown)

Planar head support

Micro adjust head support

Kid Kart swing away lateral support - push button release

TLC swing away lateral support - spring release

Adjustable contour cushion (shown)

- Centrally adjustable pre-ischial ledge
- Undercut front ledge
- Adjustable foam wedges for lateral thigh and knee supports
- Adjustable foam wedge for medial knee/thigh support
- Ability to attach lap belt at 45° or 90°
- Seat well foam options conform med, med/soft or soft

Butterfly harness

Lap belt

Optional hard mounted lateral hip and knee supports for adjustable contour cushion

Optional medial knee support

Footstraps

Standard seat cushion

- Medium density contoured foam
- Built in 3" lateral hip & thigh supports
- Medial thigh contour
- Fixed ischial ledge

All seating covers are four-way stretch, latex-free, water-resistant covers.

CRAYON

The Power Pak

Power-4-KIDS

Helps Encourage:

- Cognitive Development
- Motor Skills
- Early Independence
- Creativity
- Social skills
- Fun!



Convenient, removable push handles

Center-mount joystick accommodates left hand or right hand control

Bus tie down locations tested and approved for transit

Storage for joystick prevents damage when Xpress seating shell is removed


Utilizes Kid Kart Xpress Seating Shell and the Quickie® Z-500 Power Base

Accessories



Ventbase

Ventbase with separately mounted vent and battery trays.




Upper extremity support tray



I.V. Pole/O₂ Tank Holder Option and Vent Tray

Standard vent/battery tray.




Kid Kart Specifications




Xpress

TLC

CLIENT MEASUREMENTS		
Top of head to seat	13" - 22"	15" - 25"
Top of shoulder to seat	7.5" - 16"	10" - 20"
Under arm to seat	4.5" - 12"	9" - 17"
Trunk depth	3" - 7"	8- 11"
Seat depth	5.5" - 12"	5" - 14"
Head width	3.75" - 7"	3.75" - 7"
Shoulder width	5" - 11"	7" - 12.5"
Chest width	4" - 10"	7 - 11.5"
Hip width	6" - 11.5"	6" - 13.5"
Lower leg length	5" - 12"	6.5" - 14"
PRODUCT SPECIFICATIONS		
Tilt Range	20° anterior 45° posterior	0° to 30°
Recline Range	45°	85° - 130°
Footbed Adjustment	70° to 170°	90° to 160°
Weight	34 lbs	43 lbs
Shell only	16 lbs	22 lbs
Base only	18 lbs	21 lbs
Width	24"	26.5"
Handle bar height	36" to 41"	37" to 41"
Folded Dimensions	24" x 28" x 20"	
Shell only	16" x 19" x 14"	12.5" x 18" x 28"
Base only	24" x 28" x 15"	26.5" x 14" x 28"
Booster base height range	10.5" - 24"	N/A on TLC
Feeder base height	24"	24"





We Love Kid Kart

