





CHOICE MATTERS AND THE NEW TILITE **AERO Z** DELIVERS. TIFIT JUST FOR YOU!

TILITE INNOVATION.

· Tru-Fit System:

Reduces the number of holes in the frame, improving durability. Allows for easy adjustment and reduces set-up time.

- Adjustable Camber Mount incorporates TiLite's Tru-Fit System with up to 3.5" of Rear Seat Height adjustment.
- Speedloader Adjustable Front Caster Housing (patent pending) quickly and accurately adjusts fork angle in 1° increments.
- · SlipStream Single-Sided Fork with Integrated Bumper (patent pending).
- Performance 5-Spoke Soft Roll Wheel weighs just 1.6 ounces, 40% lighter than the original 3-Spoke Wheel.
- TiShaft back release bar makes folding and unfolding the backrest easy from any angle.

TECHNICAL SPECIFICATIONS.

Frame Type Aluminum Mono-Tube Rigid - Full Adjustability

Frame Styles SuperLite Edition
Active Edition
Rugged Edition

Armed Forces Edition
Power Adaptable Frame

Total Chair Width*

Weight Limit

Seat Widths

Seat Widths

Front Seat Widths (Seat Taper)

Seat Depths

Front Seat Heights

Rear Seat Heights

Front Frame Angles

23"

265 lbs.

3 Options

3 Options

10" - 20"

16" - 21"

12.5" - 20"

75°, 80°, 85° or 90°

Footrest Width Options Standard Taper or V-Front End O°, 2°, 4°, 6°, 8°, 12°, 15° or 18°

Front Wheels 22 Options
Rear Wheels 62 Options
Handrims 12 Options
Center of Gravity Adjustment 6"

Transport Weight* 12.4 lbs.

Convenience Transport Bracket Option

Folding Angle Adjustable Seat Back Depth Adjustable Seat Back

Aesthetics 6 Color Anodize Packages

19 Frame Colors 7 Ultrasuede Colors 8 Spinergy Wheel Colors

HCPCS Code K0005

Warranty Lifetime on Frame; Limited Lifetime on Components

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.6 lbs.

TIFIT IS EXCLUSIVELY OFFERED BY TILITE.

The difference between a TiFit chair and all other chairs is that a TiFit chair is tailored specifically to you: Every frame is made to measure. No tube is cut, bent, or welded until our skilled fabricators have a full-scale drawing of your chair based on your exact measurements. Only TiLite offers TiFit. The optimal fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION THE TILITE WAY.

Custom configuration is when you are able to select options and components on a chair to meet your needs. While other companies offer some custom configuration, TiLite offers the widest range of options from which to choose. More choices mean more ability to truly create a custom configuration of options that best meets your needs. Choice matters and that is why TiLite offers more options than any other manufacturer! The choice is yours.



Ti

TiLite | 800.545.2266 | tilite.com © 2011 TiSport, LLC.

All TiLite wheel chairs are shown with available options. CF_AEROZS2_0911RevB