







All TiLite wheel chairs are shown with available options, CF TR 0911RevC

THE TILITE TR AND YR. RAISING THE BAR TO A WHOLE NEW LEVEL. LITERALLY!

WHY TITANIUM?

Reason No. 1: Titanium Absorbs Vibration. Reason No. 2: Titanium Won't Break.

Reason No. 3: Titanium has the Highest Strength-to-Weight Ratio of Any Metal on Earth.

TILITE INNOVATION.

- · Compact Dual-Tube 1" aerospace grade titanium frame optimizes energy transfer for the most efficient propulsion.
- · Bullet caster housing reduces weight and gives the TR and YR a stylish look.
- Performance 5-Spoke Soft Roll Wheel weighs just 1.6 ounces, 40% lighter than the original 3-Spoke Wheel.
- TiShaft back release bar makes folding and unfolding the backrest easy from any angle.

TECHNICAL SPECIFICATIONS.

Titanium Dual-Tube Rigid - Minimal Adjustability Frame Type

Frame Style Power Adaptable Frame Total Chair Width* TR 23" and YR 19"

250 lbs. Weight Limit Seat Widths 12" - 20" Front Seat Widths (Seat Taper) 4 Options Seat Depths 12" - 20"

Ergonomic Seats 7 Options Front Seat Heights 15" - 21" Rear Seat Heights 12.5" - 21'

Front Frame Angles 60°, 65°, 70°, 75°, 80°, 85° or 90° Footrest Width Options Standard Taper or V-Front End 0°, 2°, 4°, 6°, 8°, 12°, 15° or 18° Cambers

Front Wheels 22 Options Rear Wheels 63 Options Handrims 12 Options Center of Gravity Adjustment

Transport Weight* TR 9.9 lbs. and YR 8.8 lbs. Folding Angle Adjustable Seat Back Convenience Aesthetics

19 Frame Colors

4 Color Anodize Packages 7 Upholstery Colors 8 Spinergy Wheel Colors

6 Titanium Frame Finishes

HCPCS Code TR (K0009) and YR (E1229) Warranty Lifetime on Frame; Limited Lifetime on Components

eights and widths based on 16" x 16" TR and 12" x 12" YR chair with 0" camber, minimum rear wheel spacing and all standard ures. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.6 lbs. and a pair of rear wheels weighs 7 lbs.

TIFIT IS EXCLUSIVELY OFFERED BY TILITE.

The difference between a TiFit chair and all other chairs is that a TiFit chair is tailored specifically to you: Every frame is made to measure. No tube is cut, bent, or welded until our skilled fabricators have a full-scale drawing of your chair based on your exact measurements. Only TiLite offers TiFit. The optimal fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION THE TILITE WAY.

Custom configuration is when you are able to select options and components on a chair to meet your needs. While other companies offer some custom configuration, TiLite offers the widest range of options from which to choose. More choices mean more ability to truly create a custom configuration of options that best meets your needs. Choice matters and that is why TiLite offers more options than any other manufacturer! The choice is yours.