

QUANTUM[®] 1450

The Quantum® 1450 front-wheel drive bariatric power base features innovative design that provides exceptional outdoor performance while delivering outstanding tight space maneuverability. The 1450 accepts a complete range of rehab seating and electronics options.



the spling to the spling of th

Specifications

DRIVE WHEELS	14" Solid
CASTER WHEELS	9" Solid
ANTI-TIP WHEELS	4" solid
MAXIMUM SPEED ⁴	Up to 4.5 mph
GROUND CLEARANCE	1.5" (anti-tip wheels)
TURNING RADIUS 5	26.5"
OVERALL LENGTH	38.25" without foot riggings
BASE WIDTH	29.25"
SEATING SIZES	TRU-Balance® 2 HD Tilt
Synergy®/Static Seat	Width: 20-32" Depth: 20-28"
Power Tilt	Width: 20-32" Depth: 20-28"
Lift & Tilt	N/A
SEAT-TO-FLOOR HEIGHTS	TRU-Balance® 2 HD Tilt
SEAT-TO-FLOOR HEIGHTS Synergy®/Static Seat 8	TRU-Balance® 2 HD Tilt 16.5-18.5"
Synergy®/Static Seat ⁸	16.5-18.5"
Synergy®/Static Seat ⁸ Power Tilt ⁸	16.5-18.5" 18.125-20.125"
Synergy®/Static Seat ⁸ Power Tilt ⁸ Lift & Tilt ⁸	16.5-18.5" 18.125-20.125" N/A
Synergy®/Static Seat ⁸ Power Tilt ⁸ Lift & Tilt ⁸ MANUAL TILT/RECLINE ¹⁰	16.5-18.5" 18.125-20.125" N/A
Synergy®/Static Seat ⁸ Power Tilt ⁸ Lift & Tilt ⁸ MANUAL TILT/RECLINE ¹⁰ BATTERY SIZE ⁶	16.5-18.5" 18.125-20.125" N/A No Group 24
Synergy®/Static Seat 8 Power Tilt 8 Lift & Tilt 8 MANUAL TILT/RECLINE 10 BATTERY SIZE 6 BATTERY WEIGHT 14	16.5-18.5" 18.125-20.125" N/A No Group 24 52 lbs. 90A VR2 6-key
Synergy®/Static Seat 8 Power Tilt 8 Lift & Tilt 8 MANUAL TILT/RECLINE 10 BATTERY SIZE 6 BATTERY WEIGHT 14 AVAILABLE ELECTRONICS 12	16.5-18.5" 18.125-20.125" N/A No Group 24 52 lbs. 90A VR2 6-key 100A Q-Logic 2 EX
Synergy®/Static Seat 8 Power Tilt 8 Lift & Tilt 8 MANUAL TILT/RECLINE 10 BATTERY SIZE 6 BATTERY WEIGHT 14 AVAILABLE ELECTRONICS 12 BATTERY CHARGER	16.5-18.5" 18.125-20.125" N/A No Group 24 52 lbs. 90A VR2 6-key 100A Q-Logic 2 EX 8A Off-board
Synergy®/Static Seat 8 Power Tilt 8 Lift & Tilt 8 MANUAL TILT/RECLINE 10 BATTERY SIZE 6 BATTERY WEIGHT 14 AVAILABLE ELECTRONICS 12 BATTERY CHARGER MOTOR PACKAGES	16.5-18.5" 18.125-20.125" N/A No Group 24 52 lbs. 90A VR2 6-key 100A Q-Logic 2 EX 8A Off-board 4-Pole 4.5 mph

- 1) Standard
- 2) Optional
- 3) Requires HD Suspension
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems.
- 7) Weight capacity is dictated by the selection of seating.
- 8) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 9) 250 lbs. weight capacity with TRU-Balance® Systems
- 10) The Synergy® Manual Tilt weight capacity is 75 lbs. and maximum dimensions are 14" W x 14" D. Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 11) Requires Group 34 batteries
- 12) All Quantum Rehab® units can be ordered with specialty controls
- Required with high-speed motors
- 14) Battery weight may vary +/- 3 lbs. based on manufacturer.

Features

- 600 lbs. weight capacity
- Maximum speed up to 4.5 mph
- High-performance, high-efficiency motor package
- Accepts a wide range of seating and positioning systems

Options

- TRU-Balance® 2 HD Power Tilt
- Manual and Power Elevating Leg Rests
- Comfort seating
- Q-Logic 2 controller

Colors







Beamsville, Ontario • 888-570-1113 • www.quantumrehab.com • www.ilevel.rehab









