

Quickie QXi AADL



EIF6, October 2018

	Mark For:				Submit	ting for:	□ Quote	 □ 0	rder
	Date:					PO#:			
	Dealer Acct #:					•			
	Dealer:						ADDITION	AL SHIPPING	INFORMATION
	Dealer Contact:					Ship To:			
	Dealer Address:					Attention:			
	Dealer City:	PV:PC:				Address:			
	Dealer Phone: () Fax: ()				Address:			
(Confirmation Email:					ip To City:		PV:	PC:
	Confirm Via:	☐ Fax ☐ Email			Ship	To Phone:	()	Fax: ()
~	Option #	Definition	Price	~	Option #		De	finition	Price
			1 F	RAME					
		Base Model					Frame Styl	е	
	196WM78	Quickie QXi (Adjustable Axle) 265 lbs	\$2,435.00		Seat Height I				
_		For 14" sizes use E1236	ΨΞ, :σσ:σσ		196FS9		Standard Frame		NC
					196FS30		Hemi Frame		NC
							Transit		NEW!
				<u> </u>	196TO38		Transit System V	VC19	\$363.00
				-					
			2. SIZE /	SEATIN	G				
				-					
_	196W5	Frame Width 14"	NO	<u> </u>	D "	1 0 10	Seat Deptl	1	
	196W5	15"	NC NC	<i>⊢ram</i> □	e Depths equ 196SD7		<i>ертпѕ.</i> 14"		NC
	196W7	16"	NC		196SD7		15"		NC NC
	196W8	17"	NC		196SD9		16"		NC
	196W9	18"	NC		196SD10		17"		NC
	196W10	19"	\$276.00		196SD11		18"		NC
	196W11	20"	\$276.00		196SD71		19"		\$276.00
			AUELGU		196SD13		20"		\$276.00
		Seat Sling	NEW!						
	196SL1	Standard Seat Sling	NC	Ļ	1 1437		AY Cushions (O		" 1
	196SL6	Omit Seat Sling	-\$20.00				n-chair, please se		
		Made-to-Order Seating		Cush	ion name or		<i>r must be filled o</i> lled Cushion	ut with desired	wiath and depth
	196QC01	Surefit	NC		196JC20		Name or Part #		Price List
	196QC02	Configurfit	NC	_	1000020	Width:		epth:	1 1100 2101
				-					
	196QC03	JayFit	NC	-		A			

Quickie QXi featuring the Jay ION™ cushion

Quickie QXi, Ultra Lightweight Adjustable Folding Wheelchair

3. CASTERS / SEAT HEIGHT

	Front Seat Height							
See	Front Seat He	ight Matrix. 1	2.5 & 13	.0 only a	vailable w/ 60	° hanger		
12.5	" & 13" only av	ailable w/ Alu	ım Angle	Adj Foo	tplate			
	196SH34	12.5"	NC		196SH14	16.0"	NC	
	196SH35	13.0"	NC		196SH15	16.5"	NC	
	196SH40	13.5"	NC		196SH16	17.0"	NC	
	196SH41	14.0"	NC		196SH17	17.5"	NC	
	196SH42	14.5"	NC		196SH18	18.0"	NC	
	196SH12	15.0"	NC		196SH20	18.5"	NC	
	196SH13	15.5"	NC		196SH21	19.0"	NC	
					196SH23	19.5"	NC	

Caster Wheels				
196CW3	4" Low Profile Polyurethane	NC		
196CW44	4" x 1.25" Semi Pneumatic	NC		
196CW5	5" Low Profile Polyurethane	NC		
196CW7	5" x 1.5" Semi Pneumatic	\$98.00		
196CW34	6" Polyurethane	NC		
196CW35	6" x 1.5" Semi Pneumatic	\$98.00		
196CW32	8" Polyurethane	NC		
196CW39	8" x 1.5" Semi Pneumatic	\$130.00		

Caster Fork				
196CF71	3.5" (standard stem bolt only)	NC		
196CF12	5.25"	NC		
196CF17	6"	NC		
196CF23	7"	NC		

Caster Stem Bolt					
196CS1	Standard	NC			
196CS5	3/4"	NC			
196CS9	1 1/2"	NC			

Front Seat Height Matrix

Front Seat Height Matrix - Standard Stem Bolt								
		3.5" FORK				5 1/4" FORK		
Caster	S	ΓD	HE	MI	STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	14.5	15.0	12.5	13.0	16.0	17.0	14.0	15.0
4 x 1.25 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 Low Poly	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 x 1.5 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
6 Poly	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
8 Poly	N/A	N/A	N/A	N/A	19.0	19.0	17.0	17.0
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.0	16.5	17.0
	6" FORK			7" FORK				
Caster	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0
4 x 1.25 Semi Pneu	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0
5 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
5 x 1.5 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
6 Poly	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0
6 x 1.5 Semi Pneu	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0
8 Poly	19.5	19.5	17.5	18.0	19.5	19.5	17.5	19.0
8 x 1.5 Semi Pneu					18.5			
Seat squeeze a	and se	at dep	th car	affec	t FSH	availa	bility	
For 3/4" and	d 1 1/2	2" Ster	n bolt	matrix	see p	age 5	5	

4. REAR WHEELS / REAR SEAT HEIGHTS

	Rear Seat Height						
See	See Rear Seat Height Matrix						
	196RH9	12"	NC NC				
	196RH10	12.5"	NC				
	196RH11	13"	NC				
	196RH12	13.5"	NC				
	196RH13	14"	NC				
	196RH14	14.5"	NC				
	196RH15	15"	NC				
	196RH16	15.5"	NC				
	196RH17	16"	NC				
	196RH18	16.5"	NC				
	196RH19	17"	NC				
	196RH26	17.5"	NC				
	196RH20	18"	NC				
	196RH27	18.5"	NC				
	196RH21	19"	NC				

	Rear Wheel Size NEW!					
See	Rear Tire & W	heel Matrix for availability	<u>-</u> -			
	196RW3	20" Mag Wheel	NC			
	196RW8	22" Mag Wheel	NC			
	196RW1	24" Mag Wheel	NC			
	196RW103	20" Mag Wheel 5 Spoke	NC			
	196RW104	22" Mag Wheel 5 Spoke	NC			
	196RW102	24" Mag Wheel 5 Spoke	NC			
	196RW90	20" Lite Spoke Wheel	\$286.00			
	196RW91	22" Lite Spoke Wheel	\$286.00			
	196RW92	24" Lite Spoke Wheel	NC			
	196RW93	25" Lite Spoke Wheel	NC			
	196RW94	26" Lite Spoke Wheel	\$113.00			
	196RW25	Omit	-\$43.00			

Rear Seat Height Matrix					
Rear Wheel Size	Std Frame	Hemi Frame			
20" Rear Wheel	13.0" - 16.5"	12.0" - 15.0"			
22" Rear Wheel	14.0" - 17.5"	12.5" - 16.0"			
24" Rear Wheel	15.0" - 18.5"	13.5" - 17.0"			
25" Rear Wheel	15.5" - 19.0"	14.0" - 17.5"			
26" Rear Wheel	16.0" - 19.5"	14.5" - 18.0"			

	Rear Wheel Tire NEW!						
See	Rear Tire & W	heel Matrix for availability					
	196RT1	Pneumatic	NC				
	196RT2	Pneumatic w/ Airless Insert	\$221.00				
	196RT3	Full Poly	NC				
	196RT39	High-Pressure Clincher Black	\$93.00				
	196RT21	Schwalbe Marathon Plus Evolutio	n \$141.00				

	Handrim NEW!					
Natu	ral Fit only ava	ailable w/	22, 24, 25 Lite Spoke. Plastic NA w/ 25	, 26"		
See	Apendix for Re	ear wheel	/ Handrim matrix			
	196HR1		Aluminum Anodized	NC		
	196HR4		Plastic	\$130.00		
	196HR22		Natural Fit - Standard Grip	\$346.00		
	196HR23		Natural Fit - Super Grip	\$389.00		
	196HR32		Natural Fit - No Thumb Grip	\$422.00		
	196HR26		The Surge - Oval w/Gription Strip	\$446.00		
	196HR27		The Surge LT - Oval w/Gription Strip	\$488.00		
	196HR28		Q-Grip Neoprene Coated	\$403.00		

Seat Squeeze Matrix

Frame Depth	Max Seat Squeeze
14"	2.0"
15" - 16"	2.5"
17" - 18"	2.5"
19" - 20"	3.0"

4. REAR WHEELS / REAR SEAT HEIGHTS (CONT)

Rear Tire & Wheel Matrix

		el	
Tire	5 Spoke Mag	Lite Spoke	Mag Wheel
Pneumatic	20,22,24	20,22,24,26	20,22,24
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26	20,22,24
Full Poly	20,22,24	20,22,24	20,22,24
High-Pressure Clincher Blk	N/A	22,24,26	N/A
Schwalbe Marathon	N/A	24,25,26	N/A
20" = 451mm 22" = 501mm	24" = 540mr	n 25" = 559mm	26" = 590mm

			One Arm Drive	NEW!		
Avai	Available only with 24" MBL Lite spoke wheels. Std Aluminum black handrim					
	196AL16		One Arm Drive Adj. Axle Plate	\$1,239.00		
	196HM1		☐ Right Hand Mount	NC		
	196HM2		☐ Left Hand Mount	NC		

	Axle	NEW!
196AX1	Quick Release Steel	NC
196AX11	Quad Quick Release Steel NA on Spoke	\$163.00
196AX45	Threaded Axle	\$134.00

	Adjustable Axle Options	MEMH
196AS4	Adjustable sleeve	NC
196AS5	Fixed Sleeve 0°	NC
196AS6	Fixed Sleeve 3°	NC

		Center of Gravity	NEW!				
Mea	Measure from front of back post to center of rear axle						
	196CC2	1"	NC				
	196CC3	1.5"	NC				
	196CC4	2"	NC				
	196CC5	2.5"	NC				
	196CC6	3"	NC				

5. WHEEL LOCKS

	Wheel Lock	NEW!
196WL1	High Mount Push	NC
196WL2	High Mount Pull	NC
196WL30	Compact	\$141.00
196WL24	Unilateral Wheel Lock Right	\$146.00
196WL23	Unilateral Wheel Lock Left	\$286.00
196WL10	Omit Wheel Lock	\$286.00
196WL12	Do Not Mount	NC

		Extension Handle			
Ext. I	Ext. handles only available with high mount push or pull wheel locks				
	196EH7	Extension Handle (Pair)	\$65.00		

6. HANGERS AND FOOTRESTS

	Front Hangers				
Foot	olate mounts v	with Extension Tube using clamping hardware			
	196H145	60° Swing In/Out Hanger	NC		
	196H139	70° Swing In/Out Hanger	NC		
	196H141	80° Swing In/Out Hanger	NC		
	196H159	Articulating Hanger NA w/ Alum Angle Adj Front Mount footplates.	\$471.00		
	196H160	Elevating Hanger - Ext Mount NA w/ Alum Angle Adj Front Mount footplates.	\$304.00		
	196H161	Elevating Hanger - Front Mount Only available w/ front mount alum angle adj footplates and 14"-22" width	\$304.00		
	196H144	Omit Hanger & Footplate	-\$16.00		

See	Lower Leg Le	ngth Matrix for range	
	196E1	Short	NC
	196E2	Medium	NC
	196E5	Long	NC
		Footplate	NEW!
Heel	loop not inclu	ded, see Footplate Option section below	
	196F1	Composite Footplate	NC
	196F271	Composite Angle Adj-(16"-22"w)	\$157.00
	196F3	Aluminum Angle Adj. Footplate	\$157.00
	196F5	Aluminum Only available w/Elevating Legrest	206
		Hanger & Footplate Options	
	196HL12	Adjustable Heel Loop	NC
	1968A02	Gel Cover For Hangers (Pair)	\$146.00

Extension Tubes

7. BACKREST / SEATING

Back Type				
	196BT10	Standard	NC	
	196BT100	Angle Adjustable Back	\$319.00	

Backrest Height				
	196B5	Med (13"-17")	NC	
	196B9	Tall (16"-20")	NC	

		Push Handles	NEW!
Integ	gral Push Hand	lles, Folding Height Adjustable & Omit is 8° for St	andard
Back	krest, Straight i	for Angle Adj. Backrest	
	196PH6	Integral Push Handles	NC
	196PH5	Fold-Down Push Handles Only available with Adjustable Upholstery. Used in combination w/Angle Adjust only avail in 16"-20" back height	\$151.00
	196PH17	Folding Push Handle with Ht. Adj. Back	\$278.00
	196PH2	Omit Push Handles	NC

7. BACKREST / SEATING (CONT)

	Back Upholstery NEWI							
	196BU1	Standard Upholstery	NC					
	196BU88	Standard Upholstery EXO	\$67.00					
	196BU70	Backrest Tension Adj Ballistic Nylon *Reference 186BU70 for billing	\$379.00					
	196BU71	\$422.00						
	196BU89	Backrest Tension Adj EXO *Reference XXXBU89 for billing	\$402.00					
	196BU14	Omit Back Upholstery	-\$20.00					
i		JAY Backs						

JAY Back (see Order forms)

Military Service Patches						
196VU1	US Marines Patch	NC				
196VU2	US Army Patch	NC				
196VU3	US Navy Patch	NC				
196VU4	US Air Force Patch	NC				
196VU12	US Coast Guard	NC				
196VU5	Purple Heart	NC				
196VU8	Purple Heart with Marines Patch	NC				
196VU9	Purple Heart with Army Patch	NC				
196VU10	Purple Heart with Navy Patch	NC				
196VU11	Purple Heart with Air Force Patch	NC				
196VU13	Purple Heart with Coast Guard	NC				

8. ARMREST / SIDEGUARDS

Price List

	Armrest - Padded Swing Away											
Some limitations when used with the Angle Adjustable Back in the rear most setting												
Please	Please circle a setting, will default to 11" if no selection is made											
П	196AA18		Pad	ded	Swing	g-Awa	ay S⊺	ΓD (9'	' to 1	3")		\$249.00
ш	1904410	Setting		9"		10"		11"		12"		13"

196JB01

Armpad - For Height Adjustable						
196AM23	Desk Length (10")	NC				
196AM24	Full Length - Full (14")	NC				

Armrest Options

		Armrest - Height Adjustable	NEW!					
Flip E	Flip Back HA Classic NA w / 14"-15" depth							
	196AA104	Single Post Ht. Adj - Std (9.75"-14.75")	\$304.00					
	196AA105	Single Post Ht Adj - Low (8.25"-12.25")	\$304.00					
	196AA65	Flip Back Ht Adj LITE (10"-15")	\$341.00					
	196AA103	Dual Post Flip Back Ht Adj	\$304.00					
	196AA159	Flip Back Ht Adj Classic (11"-15")	\$341.00					
	1964414	Omit Armrest	NC					

☐ 1968A01 Gel Cover For Swing Away Arms (pair)							
Side Guards							
196SG33	Plastic - Standard (8.5" High)	\$168.00					
196SG29	Plastic - Low (6.5" High)	\$168.00					

9. KEY OPTIONS

Anti Tips								
	196AT1	Anti Tips (Pair)	\$108.00					
_								
Pack Accessories								
	□ 196A1 Backpack - Black/Charcoal Adult							
	1961A01 Seat Pouch Black							
Positioning Belt								
	196A141	1.5" Auto Buckle w/ Clamps	\$87.00					
□ 196A804 Positioning Belt Slide Clamps \$5								

1964A01	Quickie Black Spoke Guard	\$141.00			
1964A02	\$141.00				
	Miscellaneous				
1968A03	Miscellaneous Caddy	\$81.00			

Spoke Guards

10. COLOR

NC NC			
□ 196C2 Black (Not avail as base w/Zebra) NC □ 196C141 Black Cherry NC □ 196C10 Black Opal NC □ 196C38 Blue Opal NC □ 196C36 Candy Blue NC □ 196C14 Candy Purple NC □ 196C35 Candy Red NC □ 196C142 Electric Blue NC □ 196C143 Evergreen NC □ 196C143 Evergreen NC □ 196C144 Green Apple NC □ 196C144 Green Apple NC □ 196C147 Hot Sparkle Pink NC □ 196C148 Mauve Pink NC □ 196C103 Rootbeer NC □ 196C114 Sparkle Silver NC □ 196C145 Sunrise Orange NC			
NC			
NC			

Matte Frame Colors									
196C116	Matte Black	NC							
196C151	Matte Black Cherry	NC							
196C150	Matte Electric Blue	NC							
196C149	Matte Evergreen	NC							
196C152	Matte Purple	NC							

Personalized Frame Color						
□ 196C112 Desert Camouflage						
196C43	Mossy Oak Camouflage	\$330.00				
196C42	Stars & Stripes	\$330.00				
196C108	Zebra (Choose above base color) Std with White unless above base paint is chosen NA w/ Blk, Matte Black, Blk Opal, Titanium)	\$330.00				

Additional Options Through Built for Me

To inquire about any unique modification for Quickie products, contact: Phone: (800) 263 - 3390, Fax: (800) 561-5834



Front Seat Height Matrix Cont.

Front Seat Height Matrix - 3/4" Stem									
Front S	eat H								
		3.5"	FORK			5 1/4"	FOR	K	
Caster	STD		HE	HEMI		STD		MI	
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	N/A	N/A	N/A	N/A	17.0	18.0	14.0	15.0	
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	14.5	15.5	
5 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5	
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5	
6 Poly	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0	
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0	
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	17.0	17.0	
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	19.5	19.5	16.5	17.0	
		6" F	ORK		7" FORK				
Caster	STD		HE	MI	STD HE		MI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
4 x 1.25 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
5 Low Poly	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0	
5 x 1.5 Semi Pneu	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0	
6 Poly	18.0	19.5	16.0	17.5	19.0	19.5	17.0	18.5	
6 x 1.5 Semi Pneu	18.0	19.5	16.5	17.5	19.0	19.5	17.0	18.5	
8 Poly	N/A	N/A	18.0	18.5	N/A	N/A	18.0	19.5	
8 x 1.5 Semi Pneu	19.5	19.5	17.5	18.5	N/A	N/A	18.0	19.5	
Seat squeeze a	nd se	at dep	th car	affec	t FSH	availa	bility		

Front Seat Height Matrix - 1 1/2" Stem								
3.5" FORK 5 1/4" FORK						〈		
Caster	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5
5 Low Poly	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0
6 Poly	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
		6" F	ORK			7" F	ORK	
Caster	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5
4 x 1.25 Semi Pneu	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5
5 Low Poly	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0
5 x 1.5 Semi Pneu	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0
6 Poly	18.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0
	19.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0
6 x 1.5 Semi Pneu								40.5
6 x 1.5 Semi Pneu 8 Poly	N/A	N/A	19.0	19.5	N/A	N/A	19.0	19.5
	N/A N/A	N/A N/A	19.0 18.0	19.5 19.0	N/A N/A	N/A N/A	19.0 18.5	19.5

Lower Leg Length Range Matrix

Extension Mount Hanger Leg Lengths							
Hanger	Ext Tube	Seat-to-Footrest Range					
	EXCTUDE	Composite	Comp Angle Adj	Alum Angle Adj	Alum		
con	Short	10.25" - 14.75"	10.25" - 14.75"	9.25" - 13.25"	N/A		
60° Swing In/Out	Medium	11.00" - 16.00"	11.00" - 16.00"	9.25" - 14.75"	N/A		
	Long	12.50" - 16.00"	12.50" - 16.00"	10.75" - 16.25"	N/A		
70° Swing In/Out	Short	11.00" - 15.75"	11.00" - 15.75"	10.00" - 14.50"	N/A		
	Medium	12.00" - 16.50"	12.00" - 16.50"	10.75" - 16.50"	N/A		
	Long	13.75" - 16.50"	13.75" - 16.50"	12.25" - 16.50"	N/A		
80º Swing In/Out	Short	10.50" - 15.25"	10.50" - 15.25"	10.50" - 15.25"	N/A		
	Medium	11.75" - 16.50"	11.75" - 16.50"	11.75" - 16.50"	N/A		
	Long	13.50" - 16.50"	13.50" - 16.50"	11.75" - 16.50"	N/A		

All values are ±0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights (see "Min Allowable Front Seat Heights" matrix for minimum FSH information).

Min Allowable Front Seat Heights (for Extension Mount Hangers)							
		Min Allowable FSH					
Hanger	Ext Tube	Composite	Comp Angle Adj	Alum Angle Adj	Alum		
coll	Short	14.0	13.5	12.5	N/A		
60° Swing In/Out	Medium	14.5	14.5	12.5	N/A		
	Long	16.0	16.0	13.5	N/A		
70° Swing In/Out	Short	14.0	14.0	13.0	N/A		
	Medium	15.0	15.0	13.5	N/A		
	Long	16.5	16.5	15.0	N/A		
80° Swing In/Out	Short	14.0	14.0	13.5	N/A		
	Medium	14.5	14.5	14.5	N/A		
	Long	16.5	16.5	15.5	N/A		

Rear	Wheel	& Handrim	Matrix

	Rear Wheel			
Handrim	5 Spoke Mag	Lite Spoke	Mag Wheel	
Aluminum Anodized	20,22,24	20,22,24,26	20,22,24	
Plastic Coat	20,22,24	22,24	20,22,24	
Natural Fit - Standard Grip	24	22,24,25,26	22,24	
Natural Fit - No Thumb Grip	24	22,24,25,26	22,24	
Natural Fit - Super Grip	24	22,24,25,26	22,24	
Surge - Oval w/Gription Strip	N/A	22,24,25,26	N/A	
Surge LT - Oval w/Gription Strip	N/A	22,24,25,26	N/A	
Q Grip Neoprene Coated	N/A	24,25	N/A	



ALR & ELR Hanger Leg Lengths					
Hanger	Extension		ling to trest	Min Allowable	
		Min	Max	FSH	
	Short	9 1/2	15 1/2	11 1/2	
ALR	Med	11	17 1/2	13	
	Long	13	19 1/2	15	
ELR Ext. Mount	Short	9	15	11	
	Med	9	17	11	
	Long	10	19	12	
ELR Front	No Ext (top hole)	4 1/2	8 1/2	6 1/2	
	Short	8	11	10	
Mount	Med	10	13	12	
	Long	12	15	14	

